

Australian/New Zealand Standard™

**Evaluation and routine testing in
medical imaging departments**

**Part 3.1: Acceptance tests—Imaging
performance of X-ray equipment for
radiographic and radioscopic systems
(IEC 61223-3-1:1999, MOD)**



AS/NZS 4184.3.1:2002

This Joint Australian/New Zealand Standard was prepared by Joint Technical Committee HE-008, Diagnostic Ionizing Imaging Equipment. It was approved on behalf of the Council of Standards Australia on 18 March 2002 and on behalf of the Council of Standards New Zealand on 27 February 2002. It was published on 20 May 2002.

The following interests are represented on Committee HE-008:

Australasian College of Physical Scientists and Engineers in Medicine
Australian Dental Association
Australian Institute of Radiography
Australian Radiation Protection and Nuclear Safety Agency
Australian and New Zealand Society of Nuclear Medicine
Department of Defence (Australia)
Department of Health, W.A.
Department of Human Services, Vic.
Diagnostic Imaging Association of Australasia
Medical Industry Association of Australia
Ministry of Economic Development, New Zealand
New South Wales Department of Public Works and Services
National Radiation Laboratory, New Zealand
Queensland Health
Royal Australian and New Zealand College of Radiologists
United Dental Hospital of Sydney
University of Sydney – Faculty of Dentistry

Keeping Standards up-to-date

Standards are living documents which reflect progress in science, technology and systems. To maintain their currency, all Standards are periodically reviewed, and new editions are published. Between editions, amendments may be issued. Standards may also be withdrawn. It is important that readers assure themselves they are using a current Standard, which should include any amendments which may have been published since the Standard was purchased.

Detailed information about joint Australian/New Zealand Standards can be found by visiting the Standards Australia web site at www.standards.com.au or Standards New Zealand web site at www.standards.co.nz and looking up the relevant Standard in the on-line catalogue.

Alternatively, both organizations publish an annual printed Catalogue with full details of all current Standards. For more frequent listings or notification of revisions, amendments and withdrawals, Standards Australia and Standards New Zealand offer a number of update options. For information about these services, users should contact their respective national Standards organization.

We also welcome suggestions for improvement in our Standards, and especially encourage readers to notify us immediately of any apparent inaccuracies or ambiguities. Please address your comments to the Chief Executive of either Standards Australia International or Standards New Zealand at the address shown on the back cover.

This Standard was issued in draft form for comment as DR 99318.

Australian/New Zealand Standard™

Evaluation and routine testing in medical imaging departments

Part 3.1: Acceptance tests—Imaging performance of X-ray equipment for radiographic and radioscopic systems (IEC 61223-3-1:1999, MOD)

First published as AS/NZS 4184.3.1:2002.

COPYRIGHT

© Standards Australia/Standards New Zealand

All rights are reserved. No part of this work may be reproduced or copied in any form or by any means, electronic or mechanical, including photocopying, without the written permission of the publisher.

Jointly published by Standards Australia International Ltd, GPO Box 5420, Sydney, NSW 2001 and Standards New Zealand, Private Bag 2439, Wellington 6020

ISBN 0 7337 4493 1

PREFACE

This Standard was prepared by the Joint Standards Australia/Standards New Zealand Committee HE-008, Diagnostic Ionizing Imaging Equipment. It is one of an ever-expanding series of acceptance and performance assessment Standards for use within health care facilities.

This Standard has been reproduced, with national modifications, from IEC 61223-3-1:1999, *Evaluation and routine testing in medical imaging departments*, Part 3-1: *Acceptance imaging performance of X-ray equipment for radiographic and radioscopy systems*.

Appendix ZZ lists the variations between this Standard and IEC 61223-3-1. These changes are indicated by a rule in the margin against each Clause affected.

In the text of this Standard, the following print types are used:

- (a) Requirements, compliance with which can be tested and definitions
.....in large roman type
- (b) Explanations, advice, introductions, general statements, exceptions and references
.....in smaller roman type
- (c) Headings of sub-clauses and text specifications
.....in *italic type*
- (d) Terms used throughout the Standard, which have been defined in Clause 3 or the Index of Defined Terms (see Annex A)
.....IN SMALL CAPITALS

As this publication has been reproduced from an International Standard, the following modifications apply:

- (i) Its number does not appear on each page of text and its identity is shown only on the cover and title page.
- (ii) In the source text 'this part of IEC 61223' should read 'this Australian/New Zealand Standard'.
- (iii) A full point substitutes for a comma when referring to a decimal marker.

The terms 'normative' and 'informative' have been used in this Standard to define the application of the annex or appendix to which they apply. A 'normative' annex or appendix is an integral part of a Standard, whereas an 'informative' annex or appendix is only for information and guidance.

The references to international Standards should be replaced by references to the following Australian/New Zealand Standards.

<i>Reference to International Standard</i>		<i>Australian/New Zealand Standard</i>	
IEC		AS/NZS	
60336	X-ray tube assemblies for medical diagnosis – Characteristics of focal spots	4274	X-ray tube assemblies for medical diagnosis – Characteristics of focal spots
60471	Graphical symbols for use on equipment	—	
60471-1	Part 1: Overview and application	—	
60522	Inherent filtration of an X-ray tube assembly	—	
60580	Area exposure product meter	—	
60601	Medical electrical equipment	3200	Medical electrical equipment
60601-1	Part 1: General requirements for safety	3200.1.0	Part 1.0: General requirements for safety – Parent Standard

IEC		AS/NZS	
60601-1-3	Part 1: General requirements for safety 3. Collateral standard: General requirements for radiation protection in diagnostic X-ray equipment	3200.1.3	Part 1.3: General requirements for safety Collateral Standard: Requirements for radiation protection in diagnostic X-ray equipment
60601-2-7	Part 2: Particular requirements for the safety of high-voltage generators of diagnostic X-ray generators	3200.2.7	Part 2.7: Particular requirements for safety – High voltage generators of diagnostic X-ray generators
60601-2-28	Part 2: Particular requirements for the safety of X-ray source assemblies and X-ray tube assemblies for medical diagnosis	3200.2.28	Part 2.28: Particular requirements for safety—X-ray source assemblies and X-ray tube assemblies for medical diagnosis generators
60788	Medical radiology – Terminology	—	—
60878	Graphical symbols for electrical equipment in medical practice	4334	Graphic symbols for use on medical electrical equipment
61223	Evaluation and routine testing in medical imaging departments	4184	Evaluation and routine testing in medical imaging departments
61223-1	Part 1: General aspects	4184.1	Part 1: General aspects
61267	Medical diagnostic X-ray equipment—Radiation conditions for use in the determination of characteristics	4358	Medical diagnostic X-ray equipment Radiation conditions for use in the determination of characteristics
ISO			
2092	Light metals and their alloys—Code of designation based on chemical symbols	—	—

CONTENTS

INTRODUCTION	vi
--------------------	----

Clause

1	Scope and object	1
1.1	Scope	1
1.2	Object	1
2	Normative references	2
3	Terminology	3
3.1	Degree of requirements	3
3.2	Use of terms	3
3.3	Defined terms	3
4	General aspects of ACCEPTANCE TESTS	4
4.1	General conditions to be considered in test procedures	4
4.2	Documents and data for the tests	4
4.3	Test conditions	4
4.4	Test parameters	5
4.5	Test equipment including PHANTOMS (ATTENUATION devices) and TEST DEVICES	6
4.6	Evaluating the test results	7
5	Test methods for RADIOGRAPHY EQUIPMENT	8
5.1	Visual and functional tests	8
5.2	*X-RAY TUBE VOLTAGE	8
5.3	*TOTAL FILTRATION	9
5.4	*FOCAL SPOT of the X-RAY TUBE	9
5.5	*Limitation and indication of the extent of the X-RAY BEAM	10
5.6	*Linearity and reproducibility of TRANSMISSION KERMA or RADIATION OUTPUT	12
5.7	*ATTENUATION RATIO of material between the PATIENT and the X-RAY IMAGE RECEPTOR	13
5.8	*AUTOMATIC EXPOSURE CONTROL (AEC)	13
5.9	LINE PAIR RESOLUTION for DIRECT RADIOGRAPHY	15
5.10	* AIR KERMA area product indicator	15
6	Test methods for RADIOSCOPY EQUIPMENT	15
6.1	Visual and functional tests	15
6.2	X-RAY TUBE VOLTAGE	16
6.3	TOTAL FILTRATION	16
6.4	FOCAL SPOT of the X-RAY TUBE	16
6.5	Functioning of the AUTOMATIC EXPOSURE RATE CONTROL (AERC)	16
6.6	Limitation of the extent of the X-RAY BEAM	17

Clause

6.7	ATTENUATION RATIO of material between the PATIENT and the X-RAY IMAGE RECEPTOR	18
6.8	*AIR KERMA RATE at the ENTRANCE PLANE of the X-RAY IMAGE INTENSIFIER for RADIOSCOPY	18
6.9	*Entrance AIR KERMA RATE for RADIOSCOPY with X-RAY IMAGE INTENSIFIER	19
6.10	AIR KERMA at the ENTRANCE PLANE of the X-RAY IMAGE INTENSIFIER for CINERADIOGRAPHY or other INDIRECT RADIOGRAPHY systems (excluding digital systems)	20
6.11	Entrance AIR KERMA for CINERADIOGRAPHY or other INDIRECT RADIOGRAPHY systems (excluding digital systems)	20
6.12	*LINE PAIR RESOLUTION for RADIOSCOPY with X-RAY IMAGE INTENSIFIER and for CINERADIOGRAPHY or other INDIRECT RADIOGRAPHY systems (excluding digital systems)	21
6.13	*LOW CONTRAST RESOLUTION for RADIOSCOPY with X-RAY IMAGE INTENSIFIER and for CINERADIOGRAPHY or other INDIRECT RADIOGRAPHY systems (excluding digital systems)	22
6.14	AIR KERMA area product indicator	23
7	Additional tests required for TOMOGRAPHY EQUIPMENT	23
7.1	*Requirements	23
7.2	Test method	24
8	Test report and statement of compliance	24
Annex A (normative)	Terminology – Index of defined terms	31
Annex B (informative)	Test parameters, symbols and units	34
Annex C (informative)	Examples of low contrast TEST DEVICES	35
Annex D (informative)	Examples of requirements (accuracy, tolerances, discrepancies) according to actual IEC standards or state of the art	37
Annex E (informative)	Bibliography	41
Tables		
B.1	Test parameters, symbols and units	34
D.1	Typical values of FOCAL SPOT dimensions for NOMINAL FOCAL SPOT VALUES	37
D.2	Values for the discrepancy parameters X, Y and Z according to IEC 60601-1-3	38
D.3	Typical values for the TRANSMISSION KERMA index	38
D.4	Typical values for the ATTENUATION RATIO of material between the PATIENT and the X-RAY IMAGE RECEPTOR	39
Figures		
1	Measuring arrangement for RADIOGRAPHY and RADIOSCOPY EQUIPMENT for AIR KERMA measurements	26
2	Measuring arrangement for RADIOGRAPHY and RADIOSCOPY EQUIPMENT to test geometry and resolutions	27
3	LINE PAIR RESOLUTION TEST DEVICE	28
4	TOMOGRAPHY LINE PAIR RESOLUTION TEST DEVICE	29
5	Discrepancies in visual indication of the X-RAY FIELD	30
6	Discrepancies in covering the IMAGE RECEPTION AREA	30
Appendix ZZ	Variations to IEC 61223-3-1:1999 for application in Australia and New Zealand	43

INTRODUCTION

This standard is part of a series of International Standards which give methods of acceptance testing and constancy testing for subsystems and systems (for example diagnostic X-RAY EQUIPMENT), including film processing, used in medical imaging departments.

Some provisions or statements in this standard require additional information. Such information is presented in annex D. An asterisk in the left margin of a clause or subclause indicates the presence of such additional information.

AUSTRALIAN/NEW ZEALAND STANDARD

Evaluation and routine testing in medical imaging departments

Part 3.1:

Acceptance tests—Imaging performance of X-ray equipment for radiographic and radioscopy systems (IEC 61223-3-1:1999, MOD)

1 Scope and object

1.1 Scope

This part of IEC 61223 applies to those components of X-RAY EQUIPMENT which influence the image quality and PATIENT dose of diagnostic X-ray systems using radiographic and radioscopy imaging systems.

This standard applies to the performance of X-RAY EQUIPMENT in the ACCEPTANCE TEST on the following medical diagnostic X-RAY EQUIPMENT and ASSOCIATED EQUIPMENT:

- radiography equipment, for example:
 - stationary radiography EQUIPMENT;
 - mobile radiography EQUIPMENT;
 - skull radiography EQUIPMENT;
 - lung radiography EQUIPMENT;
 - TOMOGRAPHY EQUIPMENT – excluding COMPUTED TOMOGRAPHY;
 - radiography devices (SPOTFILM DEVICES) in RADIOSCOPY EQUIPMENT;
 - angiography EQUIPMENT (excluding DSA function);
 - CINERADIOGRAPHY equipment;
- RADIOSCOPY EQUIPMENT, including:
 - combined radiographic and radioscopy EQUIPMENT.

This standard applies to the generation of X-RADIATION and ACCESSORIES of digital systems. It does not apply to any digital image acquisition or image processing parts of the above mentioned diagnostic X-RAY EQUIPMENT.

NOTE – Since the characterization of digital detectors and image processing is still under development, this will be included in a later edition of this standard.

This standard does not apply to mammographic X-RAY EQUIPMENT, RADIOTHERAPY simulators, nor to dental X-RAY EQUIPMENT.

1.2 Object

This standard defines:

- a) the parameters which describe the performance of X-RAY EQUIPMENT with regard to imaging properties and PATIENT dose;
- b) methods of testing whether measured quantities related to those parameters comply with the specified tolerances.

These methods rely mainly on non-invasive measurements, using appropriate test equipment, performed during or after the installation is completed. Signed statements covering steps of product testing at the MANUFACTURER'S site or during the installation procedure can be used as part of the acceptance testing.

The aim is to verify compliance of the installation with specifications relating to the image quality and PATIENT dose, and to detect malfunctions that are not in agreement with those specifications.

This standard does not specify tolerances for the parameters under investigation. Nor does it consider:

- c) aspects of mechanical and electrical safety,
- d) aspects of mechanical, electrical and software performance unless they are essential to the performance of the tests directly affecting image quality and PATIENT dose.

2 Normative references

The following normative documents contain provisions which, through reference in this text, constitute provisions of this International Standard. For dated references, subsequent amendments to, or revisions of, any of these publications do not apply. However, parties to agreements based on this International Standard are encouraged to investigate the possibility of applying the most recent editions of the normative documents indicated below. For undated references, the latest edition of the normative document referred to applies. Members of IEC and ISO maintain registers of currently valid International Standards.

IEC 60336:1993, *X-ray tube assemblies for medical diagnosis – Characteristics of focal spots*

IEC 60417-1:1998, *Graphical symbols for use on equipment – Part 1: Overview and application*

IEC 60522:1976, *Inherent filtration of an X-ray tube assembly*

IEC 60580:1977, *Area exposure product meter*

IEC 60601-1:1988, *Medical electrical equipment – Part 1: General requirements for safety*

IEC 60601-1-3:1994, *Medical electrical equipment – Part 1: General requirements for safety – 3. Collateral standard: General requirements for radiation protection in diagnostic X-ray equipment*

IEC 60601-2-7:1998, *Medical electrical equipment – Part 2: Particular requirements for the safety of high-voltage generators of diagnostic X-ray generators*

IEC 60601-2-28:1993, *Medical electrical equipment – Part 2: Particular requirements for the safety of X-ray source assemblies and X-ray tube assemblies for medical diagnosis*

IEC 60788:1984, *Medical radiology – Terminology*

IEC 60878:1988, *Graphical symbols for electrical equipment in medical practice*

IEC 61223-1:1993, *Evaluation and routine testing in medical imaging departments – Part 1: General aspects*

IEC 61267:1994, *Medical diagnostic X-ray equipment – Radiation conditions for use in the determination of characteristics*

ISO 2092:1981, *Light metals and their alloys – Code of designation based on chemical symbols*

3 Terminology

3.1 Degree of requirements

In this standard, certain terms which are not printed in SMALL CAPITALS have particular meanings, as follows:

- "shall" indicates a requirement that is mandatory for compliance;
- "should" indicates a strong recommendation that is not mandatory for compliance;
- "may" indicates a permitted manner of complying with a requirement or of avoiding the need to comply;
- "specific" is used to indicate definitive information stated in this standard or referenced in other standards, usually concerning particular operating conditions, test arrangements or values connected with compliance;
- "specified" is used to indicate definitive information stated by the MANUFACTURER in ACCOMPANYING DOCUMENTS or in other documentation relating to the EQUIPMENT under consideration, usually concerning its intended purposes, or the parameters or conditions associated with its use or with testing to determine compliance.

3.2 Use of terms

In this standard, terms printed in SMALL CAPITALS are used as defined in IEC 60601-1, IEC 60788, IEC 61223-1 and in 3.3 of this standard (see annex A).

NOTE – Attention is drawn to the fact, that, in cases where the concept addressed is not strongly confined to the definition given in one of the publications listed above, a corresponding term is printed in lower case letters.

3.3 Defined terms

3.3.1

ARTIFACT

apparent structure visible in the image which does not represent a structure within the object and which cannot be explained by noise or the MODULATION TRANSFER FUNCTION of the system

3.3.2

LINE PAIR RESOLUTION

highest spatial frequency of the specified line-group test pattern imaged under specified conditions which is distinguishable in the image. The unit is lp/mm

NOTE – LINE PAIR RESOLUTION is used here as a practical substitute for spatial resolution.

3.3.3

LOW CONTRAST RESOLUTION

lowest contrast detail object of a specified shape and area that can be resolved from an uniform background

3.3.4

RADIATION OUTPUT

AIR KERMA per CURRENT TIME PRODUCT (mGy/mAs) at a given distance from the FOCAL SPOT in the primary X-RAY BEAM

3.3.5

TRANSMISSION KERMA (TRANSMISSION KERMA RATE)

AIR KERMA (AIR KERMA RATE) in the central X-RAY BEAM behind the specified attenuating layer

4 General aspects of ACCEPTANCE TESTS

4.1 General conditions to be considered in test procedures

The aim of an ACCEPTANCE TEST is to demonstrate that the specified characteristics of the equipment lie within the specified tolerances. Some requirements are enforced by legislation. Other requirements and specifications may be in the order contract, in the supplier's brochure or in other standards, for example in the IEC 60601 series.

Before any ACCEPTANCE TEST according to this standard is carried out, the EQUIPMENT has to be installed and put into service according to the set-up procedure as given in the MANUFACTURER's documentation.

An inventory of the EQUIPMENT under test, the ACCOMPANYING DOCUMENTS, and the test protocols shall be compiled. Each item shall be identified by its MODEL OR TYPE REFERENCE (type number) and SERIAL NUMBER, and the entire inventory shall be compared with the order contract.

RADIOGRAPHIC CASSETTES with INTENSIFYING SCREENS, RADIOGRAPHIC FILMS and film processing are vital parts in the imaging chain. It is the responsibility of the USER to show that these components perform in an acceptable way, based upon information given by MANUFACTURERS of RADIOGRAPHIC FILMS and INTENSIFYING SCREENS, for example with respect to sensitivity, contrast and absence of ARTIFACTS.

Non-invasive measurements are preferred for ACCEPTANCE TESTS. Whenever invasive tests are part of the programme, it shall be shown that the EQUIPMENT has been restored to its pre-test condition after the test.

4.2 Documents and data for the tests

The following documentation is required:

- statements of compliance with applicable parts of IEC 60601;
- list of EQUIPMENT or EQUIPMENT parts ordered and the actual delivery list;
- performance specification as agreed upon between the purchaser and the supplier;
- results from tests performed at the MANUFACTURER's site or during installation covering items of importance to quality, such as NOMINAL FOCAL SPOT VALUE;
- INSTRUCTIONS FOR USE, including guidance for the operation of the EQUIPMENT;
- details of the actual operating conditions under which the X-RAY EQUIPMENT is to be used in medical practice and whether this results in a limitation of the scope of the tests or of the functionality of the EQUIPMENT. If certain functions are disabled, only those used need to be tested;
- guidance as to the extent and frequency of maintenance procedures;
- reports on previous tests where applicable;
- list of agreed technical modifications in the meantime between the order contract and the ACCEPTANCE TEST.

4.3 Test conditions

Different categories of tests can be identified:

- visual inspection;
- functional tests;
- system performance;
- check of the uncertainty in the values of variables.

The measuring arrangements which may be used for performing tests are illustrated in figures 1 and 2.

Figure 1 includes the measuring arrangement for the TRANSMISSION KERMA, K_T , or the TRANSMISSION KERMA RATE, $\dot{K}_T$, and the X-RAY IMAGE RECEPTOR AIR KERMA, K_B , or the X-RAY IMAGERECEPTOR AIR KERMA RATE, $\dot{K}_B$, and the test parameters derived from them.

Figure 2 includes the measuring arrangement to test geometry and resolutions.

The arrangements in figures 1 and 2 are indicative only. The test may be carried out in the vertical or the horizontal position according to the mode of operation of the EQUIPMENT. Not every component in the figures is needed in every test.

The X-RAY FIELD size shall be the minimum size required for each measurement.

The distance between the additional attenuating layer and the detector of the KERMAMETER (KERMA RATEMETER) shall be not less than 250 mm.

NOTE – For the effects of SCATTERED RADIATION, see 4.5.4.

The tests shall yield information reasonably necessary for a demonstration of performance over the full range of OPERATOR accessible variables.

All relevant data, such as the identification of the X-RAY EQUIPMENT under test, identification of the test equipment used, geometrical set-up, operating characteristics, correction factors and test results of the ASSOCIATED EQUIPMENT (for example film, screen, processing) shall be recorded with the test results. The record shall include the identification of the location, the date and the names of the persons performing the tests.

4.4 Test parameters

The following items are subject to ACCEPTANCE TESTING:

- identification of EQUIPMENT;
- check of documents;
- visual and functional tests;
- X-RAY TUBE VOLTAGE;
- CURRENT TIME PRODUCT;
- LOADING TIME;
- limitation and indication of the extent of the X-RAY BEAM;
- FOCAL SPOT;
- TOTAL FILTRATION;
- RADIATION OUTPUT;
- TRANSMISSION KERMA (TRANSMISSION KERMA RATE);
- function of the AUTOMATIC EXPOSURE CONTROL;
- ATTENUATION RATIO;
- AIR KERMA (AIR KERMA RATE);
- AIR KERMA (AIR KERMA RATE) at the ENTRANCE PLANE of the X-RAY IMAGE INTENSIFIER;
- LINE PAIR RESOLUTION;
- LOW CONTRAST RESOLUTION;
- optical density (for AUTOMATIC EXPOSURE CONTROL test).

Annex B lists symbols and units for some of the above items.

4.5 Test equipment including PHANTOMS (ATTENUATION devices) and TEST DEVICES

4.5.1 General

Measuring equipment used for ACCEPTANCE TESTS shall be certified (for example calibration according to national or international regulations, if applicable).

The uncertainty of measuring instruments shall be less than one-third of the specified tolerances for the quantities being measured.

4.5.2 High-voltage measuring instrument

The high-voltage measuring instrument shall measure the value of X-RAY TUBE VOLTAGES within the specified range. Instruments based on either direct or indirect measurements may be used.

4.5.3 KERMAMETER

The KERMAMETER (KERMA RATEMETER) shall have a range sufficient to measure the AIR KERMA (AIR KERMA RATE) within the required accuracy for the system under test and shall be calibrated for the applied beam qualities.

If legal regulations require the use of other dosimetric quantities, they may be applied.

4.5.4 PHANTOMS (ATTENUATION devices) and TEST DEVICES

PHANTOMS and TEST DEVICES may consist of attenuating layers (PHANTOM part of the object) and/or structural elements (TEST DEVICE part of the object) which can be arranged in combination or separately.

The following requirements apply:

a) External dimensions

PHANTOM dimensions shall be at least large enough to intercept the entire RADIATION BEAM for all test conditions applicable; see figures 1 and 2.

b) ATTENUATION and hardening

The attenuating layers of the PHANTOMS shall be in aluminium of at least 99,5 % purity (Al 99,5 according to ISO 2092) and a material thickness of $25 \text{ mm} \pm 0,5 \text{ mm}$; see IEC 61267.

Some, but not all, tests will need an additional homogeneous attenuating layer of about 1,5 mm copper.

A PHANTOM of low atomic number material (for example TISSUE EQUIVALENT MATERIAL) is used to test the function of the AUTOMATIC EXPOSURE CONTROL, for example 10 cm, 15 cm or 20 cm of water.

For some tests lead layers (1 mm to 2 mm thick) are needed to make lead masks or for the shielding of direct and indirect RADIATION.

c) Effects of SCATTERED RADIATION with various measuring arrangements

In all tests care shall be taken to reduce SCATTERED RADIATION to a minimum. If it is likely that SCATTERED RADIATION will significantly affect the measurement, the correction factor shall be determined and used in calculating the results.

d) Beam limiting TEST DEVICE

The beam limiting TEST DEVICE shall comprise structural elements for testing the centring, limitation and indication of the extent of the X-RAY BEAM as marking elements and a matrix with intervals of 1 cm made of RADIATION absorbing material.

These structural elements should be of such material and arranged in such a way that the function of the AUTOMATIC EXPOSURE CONTROL is unaffected.

e) **LINE PAIR RESOLUTION TEST DEVICE**

The TEST DEVICE shall comprise line-group test patterns with a lead thickness of 0,05 mm and grid groups with local frequencies of 0,6 lp/mm to 5,0 lp/mm with a gradation of less than or equal to 20 % from group to group. The outer dimensions are for example 55 mm × 65 mm; see figure 3.

f) **LOW CONTRAST RESOLUTION TEST DEVICE**

There are many devices available to measure the LOW CONTRAST RESOLUTION. If this parameter is measured, the results should be recorded together with the description of the TEST DEVICE used.

The detail diameters shall be such that their resolution is neither enhanced nor degraded by the frequency response of the X-RAY IMAGE INTENSIFIER-television system (imaging system). For examples, see annex C.

4.5.5 Lens

A magnifying lens shall be available. A 2 × to 6 × magnification is usually suitable.

4.5.6 Densitometer

The densitometer shall cover the optical density range 0 to 3,5.

4.5.7 Additional inspection and TEST DEVICES for TOMOGRAPHY X-RAY EQUIPMENT

The following TEST DEVICES are required:

- TEST DEVICE to test the layer height adjustment:
Holder for metallic slab with holes drilled at constant intervals (or the TEST DEVICE described above with additional holes; see figure 4) with an inclination of the short axis of the test pattern of 20° to 45° to the plane of the PATIENT SUPPORT. The holes in each row of the test pattern shall be at such intervals as to give an interval of 1 mm in the direction of height;
- TOMOGRAPHY movement TEST DEVICE:
Pinhole DIAPHRAGM for displaying the tomographic movement;
- TOMOGRAPHY LINE PAIR RESOLUTION TEST DEVICE:
A line-group test pattern with a lead thickness of between 0,05 mm and 0,10 mm and grid groups with local frequencies of 0,5 lp/mm to 4 lp/mm with a gradation of less than or equal to 40 % from group to group. The outer dimensions are for example 42 mm × 110 mm; see figure 4.

4.6 Evaluating the test results

Whenever specified limiting values or tolerances are not met, verify the results by making at least two additional measurements.

In the evaluation of the results concerning limit values (upper or lower), the uncertainty in the measurement shall be taken into consideration.

5 Test methods for RADIOGRAPHY EQUIPMENT

5.1 Visual and functional tests

5.1.1 Requirements

The operation and functioning of the X-RAY EQUIPMENT shall comply with what is specified.

All OPERATOR accessible controls shall be labelled with a graphical symbol, for example IEC 60417-1 or IEC 60878, and/or in plain language. The colour of indicator lamps shall comply with applicable standards, for example IEC 60601-1. The marking on the X-RAY TUBE ASSEMBLY shall comply with IEC 60601-2-28. The LIGHT FIELD shall be recognizable in the specified ambient illumination.

The INSTRUCTIONS FOR USE shall describe comprehensively how the X-RAY EQUIPMENT under test is to be operated. The function of each OPERATOR accessible control, indicator and display shall be described and all symbols shall be illustrated with their significance. Reproductions in the INSTRUCTIONS FOR USE shall be in agreement with the actual X-RAY EQUIPMENT, with respect to position, labels and symbols. The INSTRUCTIONS FOR USE shall be written in the language that is required locally or specified in the order contract.

5.1.2 Test methods

The tests are performed by visual inspection and functional check. They comprise:

- *inventory of EQUIPMENT under test;*
- *check on the presence of all documents according to 4.2;*
- *functional test of the mechanical and electrical adjustment devices (including check of the accuracy of FOCAL SPOT TO IMAGE RECEPTOR DISTANCE indicator and any detent positions, where this parameter is variable);*
- *functional test and identification of the control elements;*
- *visual inspection of the labelling of control elements;*
- *visual inspection of the markings on the X-RAY SOURCE ASSEMBLY;*
- *visual inspection of the brightness of the LIGHT FIELD;*
- *visual inspection of the INSTRUCTIONS FOR USE (see IEC 60601-1).*

5.2 *X-RAY TUBE VOLTAGE

5.2.1 Requirements

The MEASURED VALUES of the X-RAY TUBE VOLTAGE shall agree with the values indicated at the control panel within the specified tolerances.

5.2.2 Test methods

The measurements should be carried out using the arrangement in figure 1 without the attenuating layer. The tests are preferably performed using a non-invasive method.

Place the detector of the high-voltage measuring instrument in the centre of the RADIATION BEAM.

At least three standard X-RAY TUBE VOLTAGE measurements shall normally be performed, at 60 kV, 80 kV and 100 kV or at voltages close to these values, at 50 % or more of the highest available X-RAY TUBE CURRENT and a LOADING TIME of approximately 0,1 s. At 80 kV, the measurement shall be carried out also at the lowest and the highest available X-RAY TUBE CURRENT settings.

NOTE 1 – The suggested LOADING FACTORS given above are a minimum recommendation. The choice of LOADING FACTORS should be made considering the nature of the EQUIPMENT under test and its clinical application, so as to explore enough of the X-RAY TUBE VOLTAGE/LOADING TIME/X-RAY TUBE CURRENT relationship to establish the conformance of the EQUIPMENT to its specification and the needs of the USER.

NOTE 2 – At no time during testing should the maximum ratings of the X-RAY TUBE ASSEMBLY be exceeded.

Compare the MEASURED VALUES of the X-RAY TUBE VOLTAGE with the values indicated at the control panel and the specified tolerances.

NOTE – The LOADING FACTORS should be recorded because the X-RAY TUBE VOLTAGE measured may vary with the LOADING FACTOR.

5.3 *TOTAL FILTRATION

5.3.1 Requirements

The minimum TOTAL FILTRATION shall be as specified.

5.3.2 Test method

Compliance with the specification is checked by inspection of the markings on the X-RAY SOURCE ASSEMBLY and by examination of the ACCOMPANYING DOCUMENTS. The QUALITY EQUIVALENT FILTRATION may be determined according to clauses 3 and 4 of IEC 60522, if necessary.

NOTE – This requires measuring the HALF-VALUE LAYER under NARROW BEAM CONDITIONS with the X-RAY EQUIPMENT operating at appropriate values of X-RAY TUBE VOLTAGE and corresponding LOADING FACTORS, and to compare with the HALF-VALUE LAYER from an X-RAY TUBE with the same TARGET material and TARGET ANGLE.

A simplified HALF-VALUE LAYER measurement may be carried out with the arrangement in figure 1, without the attenuating layer. Measure the first HALF-VALUE LAYER with the X-RAY EQUIPMENT operating at appropriate values of X-RAY TUBE VOLTAGE with suitable LOADING FACTORS in the range of NORMAL USE. This test gives only an approximate estimate of the TOTAL FILTRATION because these test conditions do not fully comply with IEC 60522.

State either the compliance of the TOTAL FILTRATION with the specification, or the compliance of the measured HALF-VALUE LAYER with the requirements according to IEC 60601-1-3.

5.4 *FOCAL SPOT of the X-RAY TUBE

5.4.1 Requirements

The actual FOCAL SPOT dimensions for the stated NOMINAL FOCAL SPOT VALUES shall comply with the dimensions specified in IEC 60336. Additional specifications, for example concerning dimensions, direction of the REFERENCE AXIS or LOADING FACTORS, are subject to testing within the scope of this standard only if these specifications also state the test method.

5.4.2 Test method

The compliance of the ACTUAL FOCAL SPOT dimensions for the stated NOMINAL FOCAL SPOT VALUES with IEC 60336 shall be confirmed by the MANUFACTURER.

NOTE – FOCAL SPOT measuring procedures by SLIT CAMERA, PINHOLE CAMERA, star pattern evaluation and Fourier transform of images of TEST DEVICES all give different results concerning size and resolution. The standard FOCAL SPOT measurement is specified according to IEC 60336 by SLIT CAMERA under specified projection conditions and optical density.

5.5 *Limitation and indication of the extent of the X-RAY BEAM

5.5.1 Accuracy of marked and written indications of the X-RAY FIELD size

5.5.1.1 Requirements

The actual size of the X-RAY FIELD shall comply with what is indicated on the EQUIPMENT within the specified tolerances.

5.5.1.2 Test methods

Compliance is checked by inspection of the X-RAY EQUIPMENT and by examination of the ACCOMPANYING DOCUMENTS. Where appropriate, measure the dimensions of the X-RAY FIELD along its two major axes at selected indicated settings of the BEAM LIMITING SYSTEM, and the FOCAL SPOT TO IMAGE RECEPTOR DISTANCE, as available for NORMAL USE and at the commonly used X-RAY BEAM angulation.

Other test procedures than those described here may be used if they lead to comparable results.

Use the FOCAL SPOT TO IMAGE RECEPTOR DISTANCE at the value indicated on the EQUIPMENT or stated in the ACCOMPANYING DOCUMENTS, for the setting used.

Make the measurements with a RADIOGRAPHIC CASSETTE arranged in the tray or on the top of the PATIENT SUPPORT in the centre of the X-RAY FIELD and using the arrangement in figure 2, without the attenuating layer.

Produce RADIOGRAMS with two field sizes, for example 18 cm × 24 cm and 24 cm × 30 cm, under these conditions, with LOADING FACTORS so as to give an optical density D in the range 0,5 to 2,0 on the RADIOGRAPHIC FILM.

Measure the X-RAY FIELD size on the processed RADIOGRAPHIC FILMS and note the discrepancies with the displayed indication(s) on the X-RAY EQUIPMENT. If the RADIOGRAPHIC CASSETTES had been exposed on the top of the PATIENT SUPPORT, this field size must be corrected to correspond to that at the cassette tray.

Discrepancies between the measured X-RAY FIELD size and the displayed indications shall be within the specified tolerances.

5.5.2 Accuracy of indication of the LIGHT FIELD-INDICATOR

5.5.2.1 Requirements

The discrepancies in the indication of the LIGHT FIELD and the X-RAY FIELD shall comply with the tolerances specified.

5.5.2.2 Test methods

Site a RADIOGRAPHIC CASSETTE, for example 24 cm × 30 cm, on top of the PATIENT SUPPORT in the centre of the X-RAY FIELD and manually set a LIGHT FIELD of for example 18 cm × 24 cm. Mark the corners of the LIGHT FIELD with radio-opaque markers, for example metal wires. Choose LOADING FACTORS so as to give an optical density D in the range 0,5 to 2,0 on the RADIOGRAPHIC FILM.

Measure the X-RAY FIELD limits on the processed RADIOGRAPHIC FILM and the discrepancies from the indicated LIGHT FIELD.