



Designation: C1854 – 17

Standard Test Method for Determination of Hydrogen (total from all sources) in Mixed Oxide ((U, Pu)O₂) Sintered Pellets by the Inert Gas Fusion Technique Followed by Thermal Conductivity Measurement¹

This standard is issued under the fixed designation C1854; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ε) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 This test method covers the determination of hydrogen in nuclear-grade mixed oxides of uranium and plutonium ((U, Pu)O₂) sintered fuel pellets. This test method is an alternative to Test Method C698 for the determination of moisture in nuclear-grade sintered mixed oxide (MOX) fuel pellets. Test Method C698 describes the detection of moisture in mixed oxides using a coulometric, electrolytic moisture analyzer. Although the main source of H₂ in the fuel pellets is moisture, there could be other sources. The MOX pellet Specification C833 specifies a limit for hydrogen from all sources, not only moisture. The inert gas fusion followed by thermal conductivity detector specified in this test method allows for detection of hydrogen from all sources. Therefore, this test method can be used to determine the limit specified in C833.

1.2 The values stated in SI units are to be regarded as standard. No other units of measurement are included in this standard.

1.3 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.*

1.4 *This international standard was developed in accordance with internationally recognized principles on standardization established in the Decision on Principles for the Development of International Standards, Guides and Recommendations issued by the World Trade Organization Technical Barriers to Trade (TBT) Committee.*

2. Referenced Documents

2.1 *ASTM Standards:*²

C698 Test Methods for Chemical, Mass Spectrometric, and Spectrochemical Analysis of Nuclear-Grade Mixed Oxides ((U, Pu)O₂)

C753 Specification for Nuclear-Grade, Sinterable Uranium Dioxide Powder

C757 Specification for Nuclear-Grade Plutonium Dioxide Powder for Light Water Reactors

C833 Specification for Sintered (Uranium-Plutonium) Dioxide Pellets for Light Water Reactors

C859 Terminology Relating to Nuclear Materials

C1068 Guide for Qualification of Measurement Methods by a Laboratory Within the Nuclear Industry

3. Terminology

3.1 For definitions of terms used in this test method but not defined herein, refer to Terminology C859.

3.2 *Definitions of Terms Specific to This Standard:*

3.2.1 *MOX*—nuclear fuel composed of a mixture of uranium and plutonium oxides ((U, Pu)O₂).

3.2.2 *reference material*—material traceable to a reference material from a national standards body such as the U.S. National Institute for Standards and Technology (NIST) or equivalent.

3.2.3 *sintering*—to increase the bonding in a mass of powder or a compact by heating below the melting point of the main constituent.

3.3 *Acronyms:*

3.3.1 *LIMS*—Laboratory Information Management System

3.3.2 *TCD*—Thermal Conductivity Detector

4. Summary of Test Method

4.1 The method for the determination of total hydrogen (H₂) from all sources presented in this test method consists of fusion

¹ This test method is under the jurisdiction of ASTM Committee C26 on Nuclear Fuel Cycle and is the direct responsibility of Subcommittee C26.05 on Methods of Test.

Current edition approved June 1, 2017. Published June 2017. DOI: 10.1520/C1854-17.

² For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

of the sample in an electrode impulse furnace in a stream of nitrogen (N₂) or argon (Ar) gas at a temperature sufficient to release all hydrogen in the sample. The stream of gas carries the released hydrogen through a series of filters to remove interfering impurities and the hydrogen is measured by TCD.

4.2 The specimen, a single MOX pellet weighing approximately 6 g, contained in a small single-use graphite crucible, is fused under a flowing nitrogen (N₂) or argon (Ar) atmosphere and at a temperature greater than 1770°C. Hydrogen present in the sample is released as molecular hydrogen into the flowing carrier gas. Other gases are also liberated into the carrier gas, such as carbon monoxide, and need to be removed before they reach the thermal conductivity detector as they interfere with the measurement. These sample impurities are swept by the carrier gas through a series of filters to remove dust and water and then the gases flow through a quartz tube filled with Schuetze reagent (see 8.3) and activated charcoal. The Schuetze reagent acts as an oxidizing agent to convert carbon monoxide (CO) to carbon dioxide (CO₂). The gas then flows through a molecular sieve to separate nitrogen from hydrogen by size, delaying the larger nitrogen molecule and allowing the hydrogen to reach the TCD before the nitrogen.

4.3 The carrier gas, now free of interfering species then flows through the measurement branch of the thermal conductivity cell which is where the quantitative detection of the hydrogen released from the sample takes place. The hydrogen concentration in the carrier gas is measured by the TCD which detects small changes in the thermal conductivity of the carrier gas containing the liberated hydrogen gas compared to the thermal conductivity of the carrier gas alone.

4.4 The detector signal plotted versus time is a function of the concentration of the H₂ in a carrier gas. The area below the curve (integral) corresponds to the total amount of hydrogen in the sample. The peak is integrated by the software and the concentration is calculated taking into account the calibration factor, the blank analysis and the sample weight. The calibration of the analyzer is made by means of a reference material. Blank values are obtained from analyzing the empty crucibles. The blank results are stored. The final sample result is corrected by the blank value and the results are expressed in µg hydrogen/g MOX.

5. Significance and Use

5.1 MOX is used as a nuclear-reactor fuel. This test method is designed to determine whether the hydrogen content of the pellets meet the requirements of fuel specification. Examples of these requirements are given in Specification C833. Other requirements may apply based on agreements between the supplier and the customer.

5.2 This method is suitable for all sintered MOX pellets containing up to 15 weight % PuO₂ when the UO₂ and PuO₂ meet the requirements of Specifications C753 and C757. The method uncertainty is related to the concentration of the hydrogen in the sample. At lower concentrations, the relative

uncertainty increases. At hydrogen contents close to the typical hydrogen content specification limit (1.3 µg hydrogen/g U + Pu metal); the combined relative uncertainty at the 95 % confidence level (k = 2) is approximately 30 %.

6. Interferences

6.1 After the sample fusion, the carrier gas containing the carbon dioxide and other potentially interfering impurities (sulfur, water, and small particulate matter) pass through a series of filters and purifying reagents that remove these impurities from the carrier gas stream leaving only H₂ in the stream at the detector. If the MOX pellets are made from UO₂ and PuO₂ that meet the requirements of Specifications C753 and C757, all interferences are eliminated by the purification system.

6.2 The crucibles, if they contain hydrogen, will yield erroneously high results for the sample. The analytical method requires running blanks to correct for this potential interference in the calculation of the results (see 12.2).

6.3 The nitrogen (N₂) or argon (Ar) carrier gas could contain water and CO₂ and is filtered prior to injection in the sample combustion chamber to remove these potentially interfering components.

6.4 Weighing uncertainty of the samples is critical to the method. If the balance meets the specification in 7.1, is calibrated in accordance with manufacturer's guidance, and is checked by procedure, the potential for the balance to be a source of error is insignificant.

6.5 When using nitrogen gas with graphite crucible a secondary reaction occurs at high temperatures (>3000°C) where HCN is created, causing a bias that increases with temperature. Argon as carrier gas is used at higher temperatures to avoid this reaction.

7. Apparatus

7.1 *Analytical Balance*, with precision ±0.1 mg.

7.2 *Graphite Crucibles*—Use the crucibles of size recommended by the manufacturer of the instrument. Crucibles shall be composed of high purity graphite.

7.3 *Hydrogen analyzer*, consisting of an electrode impulse furnace suitable for operation at 1770 to 2200°C, a TCD for measuring H₂, a nitrogen or argon carrier gas injection system, auxiliary gas purification systems, and a chiller to cool the electrodes. The analyzer typically has an integrated data collection and analysis software system that allows for efficient collection and analysis of data.

7.4 *Tongs and Forceps*, for handling crucibles.

7.5 *Stainless Steel Scoops and Spatulas*, for handling pellets and reference materials.

8. Reagents and Materials

8.1 *Purity of Reagents*—Reagent grade chemicals shall be used in all tests. Unless otherwise indicated, it is intended that