



Designation: E 2154 – 01

Standard Practice for Separation and Concentration of Ignitable Liquid Residues from Fire Debris Samples by Passive Headspace Concentration with Solid Phase Microextraction (SPME)¹

This standard is issued under the fixed designation E 2154; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ε) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 This practice describes the procedure for removing small quantities of ignitable liquid residues from samples of fire debris. An adsorbent material is used to extract the residue from the static headspace above the sample. Then, analytes are thermally desorbed in the injection port of the gas chromatograph (GC).

1.2 This practice is best suited for screening fire debris samples to assess relative ignitable liquid concentration and for extracting ignitable liquid from aqueous samples.

1.3 This practice is suitable for extracting ignitable liquid residues when a high level of sensitivity is required due to a very low concentration of ignitable liquid residues in the sample.

1.3.1 Unlike other methods of separation and concentration, this method recovers a minimal amount of the ignitable residues present in the evidence, leaving residues that are suitable for subsequent resampling.

1.4 Alternate separation and concentration procedures are listed in Section 2.

1.5 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.*

2. Referenced Documents

2.1 ASTM Standards:

E 1385 Practice for Separation and Concentration of Flammable or Combustible Liquid Residues from Fire Debris Samples by Steam Distillation²

E 1386 Practice for Separation and Concentration of Flammable or Combustible Liquid Residues from Fire Debris Samples by Solvent Extraction²

E 1387 Test Method for Flammable or Combustible Liquid Residues in Extracts from Samples of Fire Debris by Gas Chromatography²

E 1388 Practice for Sampling Headspace Vapors from Fire Debris Samples²

E 1412 Practice for Separation and Concentration of Flammable or Combustible Liquid Residues from Fire Debris Samples by Passive Headspace Concentration with Activated Charcoal²

E 1413 Practice for Separation and Concentration of Flammable or Combustible Liquid Residues from Fire Debris Samples by Dynamic Headspace Concentration²

E 1459 Guide for Physical Evidence Labeling and Related Documentation²

E 1492 Practice for Receiving, Documenting, Storing, and Retrieving Evidence in a Forensic Science Laboratory²

E 1618 Guide for Ignitable Liquid Residues in Extracts from Fire Debris Samples by Gas Chromatography-Mass Spectrometry²

3. Summary of Practice

3.1 A fiber coated with a polydimethylsiloxane stationary phase is exposed to the headspace of the fire debris sample container to extract ignitable liquid residues. The fiber, which is housed in a needle similar to a syringe needle, is introduced directly in the injection port of a gas chromatograph to thermally desorb the analytes.

4. Significance and Use

4.1 This practice is suited ideally for screening samples for the presence, relative concentration, and potential class of ignitable liquid residues in fire debris.

4.2 This is a very sensitive separation procedure, capable of isolating small quantities of ignitable liquid residues from a sample, that is, a 0.1 µL spike of gasoline on a cellulose wipe inside of a 1-gal can is detectable.

4.3 Actual recovery will vary, depending on several factors, including adsorption temperature, container size, competition from the sample matrix, ignitable liquid class and relative ignitable liquid concentration.

¹ This practice is under the jurisdiction of ASTM Committee E30 on Forensic Sciences and is the direct responsibility of Subcommittee E30.01 on Criminalistics. Current edition approved Sept 10, 2001. Published December 2001.

² *Annual Book of ASTM Standards*, Vol 14.02.

4.4 Because this separation takes place in a closed container, the sample remains in approximately the same condition in which it was submitted. Repeat and interlaboratory analyses, therefore, may be possible. Since the extraction is nonexhaustive, the technique permits reanalysis of samples.

4.5 This practice is intended for use in conjunction with other extraction techniques described in Practices E 1385, E 1386, E 1388, E 1412, and E 1413.

4.6 The extract is consumed in the analysis. If a more permanent extract is desired, one of the separation practices described in Practices E 1385, E 1386, E 1412, or E 1413 should be used.

5. Apparatus

5.1 *Heating System*, such as, an oven or heating mantle to fit the evidence container (or a hot plate).

5.2 *Temperature Measuring Device*, such as, a thermometer or thermocouple capable of measuring temperatures in the range of 40 to 110°C.

5.3 *SPME Apparatus*, such as, a solid phase microextraction fiber holder.

5.4 *SPME Fiber*—Several polymer coatings can be used to extract volatiles such as ignitable liquid residues. A 100 μm polydimethylsiloxane (PDMS) coating has been shown to perform well for most C_{10} - C_{25} compounds while a 85 μm Polyacrylate (PA) and a combined 75 μm Carboxen/PDMS has been shown to perform well for C_1 - C_{10} compounds (5).

5.5 *Punch*.

5.6 *Rubber Sleeve Septum*.

6. Sample Preparation/Analysis Procedure

6.1 Observe the appropriate procedures for handling and documentation of all submitted samples as described in Practice E 1492.

6.1.1 Open and examine the fire debris sample to determine that it is consistent with its description.

6.1.1.1 Resolve any discrepancies between the submitting agent's description of the evidence and the analyst's observation with the submitting agent prior to the completion of the report.

6.2 Verify that the fiber is free of ignitable residue components and other contaminants by introducing the fiber into the GC under the injector desorption and chromatographic conditions for ignitable residue analysis. This fiber residue "blank" always should be evaluated before residue sampling.

6.3 If the sample container has a metal lid, make a hole in the lid and install a rubber sleeve septum in the opening.

6.4 Heat the container until the sample reaches the desired temperature (approximately 20-30 min, but longer equilibration times may be needed depending on the heating method or the amount and type of matrix present in the can, or both).

6.4.1 Temperatures lower than 60°C may be insufficient to volatilize compounds above C_{16} .

6.4.2 Temperatures above 80°C may result in a significant discrimination against high volatility compounds when the 100 μm PDMS fiber is used.

6.5 Remove the container from the heating apparatus.

6.6 Immediately puncture the rubber sleeve septum or plastic evidence bag with the needle of the SPME apparatus.

6.7 Expose the SPME fiber to the headspace of the container for the desired sampling duration.

6.7.1 The optimum exposure time for maximum sensitivity will depend on the temperature and the concentration and composition of the volatile compound present in the sample headspace.

6.7.2 Exposure times for routine screening of samples typically are in the range of 5-15 min.

6.7.3 Residue extracts that provide off-scale or poorly resolved chromatographic responses should be resampled at lower temperatures or shorter exposure times.

6.8 Secure the fiber in the needle by withdrawing it into the SPME apparatus. Remove the apparatus from the sample container.

6.8.1 Rubber sleeve septa are self-sealing. Puncture holes in evidence bags should be sealed after sampling.

6.9 Desorb the fiber for approximately 1.5 to 4 min at a minimum of 200°C and 260°C and analyze the adsorbate by exposing the fiber in the injection port of a gas chromatograph. Analysis is performed as described in Test Method E 1387 or Guide E 1618.

6.10 A blank shall be run before each sample to ensure that the fiber is completely desorbed and free of contaminants.

6.10.1 A blank is run by exposing and analyzing the fiber in the gas chromatograph under normal operating conditions.

6.11 Record the adsorption parameters including the exposure time and temperature in the analytical case notes.

7. Blanks and Standards

7.1 Analysis of a fiber blank (6.2) is required immediately before each sample extraction.

7.1.1 If an internal standard is routinely used, include the internal standard in the blank run by placing the internal standard in an empty container and exposing the SPME fiber as described in Section 6.

7.2 Periodically check the adsorption efficiency of the fiber by running this procedure on a sample containing a known volume of an ignitable liquid standard.

7.3 An ignitable liquid standards library should be maintained with neat and evaporated commonly ignitable liquids recovered from sample containers by this technique at various concentrations.

8. Keywords

8.1 fire debris samples; passive headspace concentration; solid phase microextraction (SPME)