

Status: Currently Official on 17-Feb-2025
Official Date: Official as of 01-Jul-2021
Document Type: USP Monographs
DocId: GUID-A0A19A5E-EFD4-42B6-B525-53CC930D70F2_4_en-US
DOI: https://doi.org/10.31003/USPNF_M4580_04_01
DOI Ref: nb00v

© 2025 USPC
Do not distribute

Benazepril Hydrochloride and Hydrochlorothiazide Tablets

DEFINITION

Benazepril Hydrochloride and Hydrochlorothiazide Tablets contain NLT 90.0% and NMT 110.0% of the labeled amount of benazepril hydrochloride ($C_{24}H_{28}N_2O_5 \cdot HCl$) and hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$).

IDENTIFICATION

- **A.** The retention times of the major peaks of the *Sample solution* correspond to those of the *Standard solution*, as obtained in the Assay.
- **B.** The UV spectra of the major peaks of the *Sample solution* correspond to those of the *Standard solution*, as obtained in the Assay.

ASSAY

• **PROCEDURE**

Buffer: Dissolve 3.45 g of [sodium phosphate, monobasic](#) in 1000 mL of [water](#). Add 2.0 mL of [triethylamine](#). Adjust with [phosphoric acid](#) to a pH of 6.8.

Mobile phase: [Acetonitrile](#) and *Buffer* (32:68)

Diluent: [Acetonitrile](#) and [water](#) (50:50)

Standard solution: Concentrations of [USP Benazepril Hydrochloride RS](#) and [USP Hydrochlorothiazide RS](#) in *Diluent*, prepared as directed in [Table 1](#). Sonicate to dissolve if needed.

Table 1

Tablet Strength Benazepril Hydrochloride/ Hydrochlorothiazide (mg/mg)	Concentration of USP Benazepril Hydrochloride RS (mg/mL)	Concentration of USP Hydrochlorothiazide RS (mg/mL)
5/6.25	0.10	0.125
10/12.5	0.10	0.125
20/25	0.10	0.125
20/12.5	0.10	0.0625

Sample stock solution: Transfer a quantity of Tablets (NLT 10) to a 200-mL volumetric flask. Add 10% of the flask volume of [water](#) and shake mechanically until the Tablets are disintegrated and well dispersed. Add another 10% of the flask volume of [acetonitrile](#), sonicate for 20 min, then mechanically shake for an additional 20 min. Dilute with *Diluent* to volume.

Sample solution: Dilute the *Sample stock solution* with *Diluent* to obtain nominal concentrations of benazepril hydrochloride and hydrochlorothiazide as directed in [Table 1](#). Pass a portion through a suitable filter and discard the first few milliliters of the filtrate.

Chromatographic system

(See [Chromatography \(621\)](#), [System Suitability](#).)

Mode: LC

Detector: UV 254 nm. For *Identification B*, use a diode array detector in the range of 210–400 nm.

Column: 4.6-mm × 15-cm; 5-μm packing [L7](#)

Flow rate: 1 mL/min

Injection volume: 10 μL

Run time: NLT 1.5 times the retention time of benazepril

System suitability

Sample: *Standard solution*

[NOTE—The relative retention times for hydrochlorothiazide and benazepril are about 0.4 and 1.0, respectively.]

Suitability requirements

Tailing factor: NMT 2.0 for benazepril and hydrochlorothiazide

Relative standard deviation: NMT 2.0% for benazepril and hydrochlorothiazide

Analysis

Samples: *Standard solution and Sample solution*

Calculate the percentage of the labeled amount of benazepril hydrochloride ($C_{24}H_{28}N_2O_5 \cdot HCl$) in the portion of Tablets taken:

$$\text{Result} = (r_U/r_S) \times (C_S/C_U) \times 100$$

r_U = peak response of benazepril from the *Sample solution*

r_S = peak response of benazepril from the *Standard solution*

C_S = concentration of [USP Benazepril Hydrochloride RS](#) in the *Standard solution* (mg/mL)

C_U = nominal concentration of benazepril hydrochloride in the *Sample solution* (mg/mL)

Calculate the percentage of the labeled amount of hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$) in the portion of Tablets taken:

$$\text{Result} = (r_U/r_S) \times (C_S/C_U) \times 100$$

r_U = peak response of hydrochlorothiazide from the *Sample solution*

r_S = peak response of hydrochlorothiazide from the *Standard solution*

C_S = concentration of [USP Hydrochlorothiazide RS](#) in the *Standard solution* (mg/mL)

C_U = nominal concentration of hydrochlorothiazide in the *Sample solution* (mg/mL)

Acceptance criteria: 90.0%–110.0%

PERFORMANCE TESTS• [DISSOLUTION \(711\)](#)**Test 1**

Medium: [0.1 N hydrochloric acid](#); 500 mL

Apparatus 1: 100 rpm

Time: 30 min

Buffer and Mobile phase: Prepare as directed in the Assay.

Standard stock solution A: 0.1 mg/mL of [USP Benazepril Hydrochloride RS](#) in *Medium*. Sonicate to dissolve if needed.

Standard stock solution B: 0.125 mg/mL of [USP Hydrochlorothiazide RS](#) in *Medium*. Sonicate to dissolve if needed.

Standard solution: [USP Benazepril Hydrochloride RS](#) and [USP Hydrochlorothiazide RS](#) in *Medium*, prepared as directed in [Table 2](#) from *Standard stock solution A* and *Standard stock solution B*

Table 2

Tablet Strength Benazepril Hydrochloride/ Hydrochlorothiazide (mg/mg)	Concentration of USP Benazepril Hydrochloride RS (mg/mL)	Concentration of USP Hydrochlorothiazide RS (mg/mL)
5/6.25	0.01	0.0125
10/12.5	0.02	0.025
20/12.5	0.04	0.025
20/25	0.04	0.05

Sample solution: Pass a portion of the solution under test through a suitable filter of 0.45-µm pore size. Discard the first few milliliters of the filtrate.

Chromatographic system

(See [Chromatography \(621\)](#), [System Suitability](#).)

Mode: LC

Detector: UV 234 nm

Column: 4.6-mm × 15-cm; 5-µm packing [L7](#)

Flow rate: 0.9 mL/min

Injection volume

For Tablet strength of 5/6.25 mg: 100 µL

For Tablet strengths of 10/12.5, 20/12.5, and 20/25 mg: 25 µL

Run time: NLT 2 times the retention time of benazepril

System suitability**Sample:** *Standard solution*

[NOTE—The relative retention times for hydrochlorothiazide and benazepril are about 0.5 and 1.0, respectively.]

Suitability requirements**Tailing factor:** NMT 2.0 for benazepril and hydrochlorothiazide**Relative standard deviation:** NMT 2.0% for benazepril and hydrochlorothiazide**Analysis****Samples:** *Standard solution and Sample solution*Calculate the percentage of the labeled amount of benazepril hydrochloride ($C_{24}H_{28}N_2O_5 \cdot HCl$) dissolved:

$$\text{Result} = (r_U/r_S) \times (C_S/L) \times V \times 100$$

 r_U = peak response of benazepril from the *Sample solution* r_S = peak response of benazepril from the *Standard solution* C_S = concentration of [USP Benazepril Hydrochloride RS](#) in the *Standard solution* (mg/mL) L = label claim (mg/Tablet) V = volume of *Medium*, 500 mLCalculate the percentage of the labeled amount of hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$) dissolved:

$$\text{Result} = (r_U/r_S) \times (C_S/L) \times V \times 100$$

 r_U = peak response of hydrochlorothiazide from the *Sample solution* r_S = peak response of hydrochlorothiazide from the *Standard solution* C_S = concentration of [USP Hydrochlorothiazide RS](#) in the *Standard solution* (mg/mL) L = label claim (mg/Tablet) V = volume of *Medium*, 500 mL**Tolerances:** NLT 80% (Q) of the labeled amount of benazepril hydrochloride ($C_{24}H_{28}N_2O_5 \cdot HCl$) and hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$) is dissolved.**Test 2:** If the product complies with this test, the labeling indicates that it meets USP *Dissolution Test 2*.**Medium:** [0.1 N hydrochloric acid](#); 500 mL**Apparatus 1:** 100 rpm**Time:** 30 min**Buffer:** 3.75 g/L of [potassium phosphate, dibasic](#) in [water](#). Adjust with [phosphoric acid](#) to a pH of 2.1.**Mobile Phase:** [Acetonitrile](#) and *Buffer* (35:65)**Standard stock solution:** 0.25 mg/mL of [USP Hydrochlorothiazide RS](#), prepared as follows. Transfer a suitable amount of [USP Hydrochlorothiazide RS](#) to a volumetric flask. Add a volume of [acetonitrile](#) equivalent to 5% of the total flask volume and, if necessary, sonicate to dissolve. Dilute with *Medium* to volume.**Standard solution:** ($L_2/500$) mg/mL of [USP Hydrochlorothiazide RS](#) from *Standard stock solution* and ($L_1/500$) mg/mL of [USP Benazepril Hydrochloride RS](#) diluted in *Medium*, where L_1 is the label claim for benazepril hydrochloride in mg/Tablet and L_2 is the label claim for hydrochlorothiazide in mg/Tablet**Sample solution:** Pass a portion of the solution under test through a suitable filter of 0.45- μ m pore size. Discard the first few milliliters of the filtrate. [NOTE—It is recommended to avoid using stainless steel cannula when sampling.]**Chromatographic system**(See [Chromatography \(621\)](#), [System Suitability](#).)**Mode:** LC**Detector:** UV 240 nm for benazepril and UV 275 nm for hydrochlorothiazide**Column:** 4.6-mm $\times$ 25-cm; 5- μ m packing [L1](#)**Flow rate:** 1 mL/min**Injection volume:** 20 μ L**Run time:** NLT 1.5 times the retention time of benazepril**System suitability****Sample:** *Standard solution***Suitability requirements****Tailing factor:** NMT 2.0 for benazepril at 240 nm and hydrochlorothiazide at 275 nm**Relative standard deviation:** NMT 2.0% for benazepril at 240 nm and hydrochlorothiazide at 275 nm