

Rubber, raw synthetic — Determination of the molecular-mass distribution of solution polymers by gel permeation chromatography

ICS 83.040.10

National foreword

This British Standard is the UK implementation of ISO 11344:2004, incorporating corrigendum April 2008.

The start and finish of text introduced or altered by corrigendum is indicated in the text by tags. Text altered by ISO corrigendum April 2008 is indicated in the text by AC1 AC1.

The UK participation in its preparation was entrusted to Technical Committee PRI/23, Test methods for rubber and non-black compounding ingredients.

A list of organizations represented on this committee can be obtained on request to its secretary.

This publication does not purport to include all the necessary provisions of a contract. Users are responsible for its correct application.

Compliance with a British Standard cannot confer immunity from legal obligations.

This British Standard was published under the authority of the Standards Policy and Strategy Committee on 31 March 2004

© BSI 2008

ISBN 978 0 580 62772 9

Amendments/corrigenda issued since publication

Date	Comments
30 September 2008	Implementation of ISO corrigendum April 2008

INTERNATIONAL STANDARD

ISO
11344

First edition
2004-04-01

Rubber, raw synthetic — Determination of the molecular-mass distribution of solution polymers by gel permeation chromatography

*Caoutchouc synthétique brut — Détermination de la répartition de la
masse moléculaire pour les caoutchoucs polymérisés en solution par
chromatographie par perméation de gel*



Contents

Page

Foreword	iv
1 Scope	1
2 Normative references	1
3 Principle	1
4 General	1
5 Reagents and materials	2
6 Apparatus	3
7 Analytical conditions	4
8 Procedure	4
9 Expression of results	7
10 Precision (only for instrumental software procedure)	8
11 Test report	9
Annex A (informative) Molecular-mass parameters determined by instrumental software	11
Annex B (informative) Calculation of molecular-mass parameters by manual procedure	15
Annex C (informative) Comparison of results obtained by automatic procedure (software) and manual procedure	19
Annex D (informative) Guidance for using precision results	20

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO 11344 was prepared by Technical Committee ISO/TC 45, *Rubber and rubber products*, Subcommittee SC 2, *Testing and analyses*.

Rubber, raw synthetic — Determination of the molecular-mass distribution of solution polymers by gel permeation chromatography

WARNING — Persons using this International Standard should be familiar with normal laboratory practice. This standard does not purport to address all of the safety problems, if any, associated with its use. It is the responsibility of the user to establish appropriate safety and health practices and to ensure compliance with any national regulatory conditions.

1 Scope

This International Standard describes a method for the determination of the molecular mass, expressed as polystyrene, and the molecular-mass distribution of polymers produced in solution which are completely soluble in tetrahydrofuran (THF) and which have a molecular-mass range from 5×10^3 to 1×10^6 .

It is not the purpose of this International Standard to explain the theory of gel permeation chromatography.

2 Normative references

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO/TR 9272, *Rubber and rubber products — Determination of precision for test method standards*

3 Principle

The molecular components of a polymer are separated on the basis of macromolecule size on a gel permeation column. A known quantity of a dilute solution of the polymer is injected into a stream of solvent, which carries it through the column at a constant rate. The concentration of the separated molecular components in the solvent stream is measured by a suitable detector. Through the use of a calibration curve, both the number-average molecular mass (M_n) and mass-average molecular mass (M_w) of the material analysed can be determined from the retention time and the corresponding concentration.

4 General

4.1 Gel permeation chromatography (GPC), which is also known as size exclusion chromatography (SEC), is a particular type of liquid chromatography which allows the separation of the various components of a polymer based on macromolecule size.

4.2 The molecules of a polymer do not all have the same mass, but comprise a range of different masses. For this reason, the usual concept of molecular mass is not applicable to polymeric materials. Instead, different average molecular masses are determined as shown in Table 1.

Table 1 — Definitions of various kinds of molecular mass

Mass-average molecular mass M_w	$\overline{[AC1]} = \Sigma(N_i M_i^2) / \Sigma(N_i M_i) \overline{[AC1]}$ $= \Sigma(A_i M_i) / \Sigma A_i$
Number-average molecular mass M_n	$= \Sigma(M_i N_i) / \Sigma N_i$ $= \Sigma A_i / \Sigma(A_i / M_i)$
z-Average molecular mass M_z	$= \Sigma(N_i M_i^3) / \Sigma(N_i M_i^2)$ $\overline{[AC1]} = \Sigma(A_i M_i^2) / \Sigma(A_i M_i) \overline{[AC1]}$
Peak molecular mass M_p	Molecular mass at peak maximum
where N_i is the number of moles having a molecular mass of M_i ; A_i is the area of the time-slice that corresponds to molecular mass M_i .	

The molecular-mass distribution is an important parameter in determining the properties of the polymer. It may be represented by the polydispersity D given by:

$$D = M_w / M_n$$

NOTE Polymers invariably consist of macromolecules with a range of molecular sizes. Even the so-called monodisperse polystyrenes have a polydispersity of 1,1 compared to a value of 1,0 for a pure compound with a single molecular mass. As the range of molecular sizes present within the polymer increases, so does the polydispersity.

5 Reagents and materials

During the analysis, unless otherwise stated, use only reagents of recognized analytical grade ~~Text~~ ~~AC1~~.

5.1 Tetrahydrofuran (THF), high-purity-grade solvent.

NOTE A large stock of THF is needed to avoid frequent refills. Changes in the quantity of dissolved air or impurities due to addition of fresh solvent cause significant variations in the refractive index and could also affect the retention time. Air bubbles at the pump head reduce the quantity of solvent pumped (leading to errors in retention volumes and times) and can block the pump if the volume of the air bubbles reaches excessive levels. After adding fresh solvent, it takes 2 to 3 hours to obtain a stable baseline.

5.2 Solution of o-dichlorobenzene (internal retention time standard) in THF, obtained by dilution of 250 mm³ (250 µl) of o-dichlorobenzene with 1 l of THF.

5.3 Set of certified polystyrene reference standards (minimum 10), with molecular masses in the range 5×10^2 to 1×10^7 (depending on the sample molecular-mass range) and a very narrow molecular-mass distribution ($D < 1,10$) (see Table 2 for an example of such a set, available from various chemical suppliers).

Table 2 — Set of polystyrene standards

Standard No.	Actual molecular mass M_i	$D (= M_w/M_n)$
1	1 030 000	1,05
2	770 000	1,04
3	336 000	1,03
4	210 000	1,03
5	156 000	1,03
6	66 000	1,03
7	30 300	1,03
8	22 000	1,03
9	11 600	1,03
10	7 000	1,04
11	5 050	1,05

6 Apparatus

Ordinary laboratory apparatus, plus the following:

6.1 Gel permeation chromatograph, consisting of the components specified in 6.1.1 to 6.1.8.

6.1.1 Solvent reservoir, of sufficient capacity to complete the analysis (see Note to 5.1) without refilling.

6.1.2 Automatic on-line degassing system or helium sparging of solvent reservoir, to stabilize the solvent flow, mainly to prevent formation of bubbles in the solvent.

6.1.3 Pump, to ensure that the THF solvent flows at a constant rate, programmable over the range 1,7 mm³/s to 165,0 mm³/s with a high degree of precision.

6.1.4 Injector or automatic sampler, with a 100 mm³ (100 µl) injection loop.

6.1.5 Columns, packed with regular, rigid, porous spheres. The pore size of column packing material is expressed in Angstrom units (1 Å = 10⁻¹⁰ m). The packing spheres are made of cross-linked polystyrene, obtained by polymerization of styrene with divinylbenzene. The spheres shall have a nominal diameter in the range 5 µm to 10 µm. The columns are generally 300 mm long. The pore size is selected depending on the range of molecular masses to be analysed.

NOTE Four columns with pore sizes 10³ Å, 10⁴ Å, 10⁴ Å and 10⁵ Å were used when the repeatability and reproducibility of the method described in this International Standard were determined. The solvent first enters the column with the lowest porosity and exits from the column with the highest porosity. Other suitable columns may be used. These types of column are available from many suppliers.

The recommended column characteristics are:

- linear range: 200 to 2 000 000;
- guaranteed column efficiency: > 50 000 plates/m;
- column arrangement: four columns (300 mm long and 4,6 mm to 8,0 mm ID).