

Whole milk — Determination of milkfat, protein and lactose content — Guidance on the operation of mid-infrared instruments

ICS 67.100 10

NO COPYING WITHOUT BSI PERMISSION EXCEPT AS PERMITTED BY COPYRIGHT LAW

BSI
Documents

AVA

Document Sharing Hub

National foreword

This British Standard reproduces verbatim ISO 9622:1999 and implements it as the UK national standard.

The UK participation in its preparation was entrusted to Technical Committee AW/5, Chemical analysis of milk and milk products, which has the responsibility to:

- aid enquirers to understand the text;
- present to the responsible international/European committee any enquiries on the interpretation, or proposals for change, and keep the UK interests informed;
- monitor related international and European developments and promulgate them in the UK.

A list of organizations represented on this committee can be obtained on request to its secretary.

Cross-references

The British Standards which implement international publications referred to in this document may be found in the BSI Standards Catalogue under the section entitled “International Standards Correspondence Index”, or by using the “Find” facility of the BSI Standards Electronic Catalogue.

A British Standard does not purport to include all the necessary provisions of a contract. Users of British Standards are responsible for their correct application.

Compliance with a British Standard does not of itself confer immunity from legal obligations.

Summary of pages

This document comprises a front cover, an inside front cover, the ISO title page, pages ii to iv, pages 1 to 27 and a back cover.

The BSI copyright notice displayed in this document indicates when the document was last issued.

This British Standard, having been prepared under the direction of the Consumer Products and Services Sector Committee, was published under the authority of the Standards Committee and comes into effect on 15 April 2000

© BSI 04-2000

ISBN 0 580 34625 0

Amendments issued since publication

Amd. No.	Date	Comments

**Whole milk — Determination of milkfat,
protein and lactose content — Guidance
on the operation of mid-infrared instruments**

*Lait entier — Détermination des teneurs en matière grasse laitière,
en protéines et en lactose — Lignes directrices pour l'utilisation
des appareils de dosage par absorption dans le moyen infrarouge*



Contents

1 Scope 1

2 Normative references 1

3 Terms and definitions 2

4 Principle..... 2

5 Principal characteristics of infrared instruments..... 2

6 Factors affecting the accuracy of measurement 3

6.1 Instrument factors 3

6.2 Physicochemical and biological factors 5

7 Calibration of the instrument..... 7

7.1 Objective..... 7

7.2 Checking the initial calibration for fat, protein and lactose 7

7.3 Maintaining calibration and confirmation of calibration validity 8

8 Sampling..... 8

9 Uniformity of test samples..... 8

10 Determination..... 8

11 Checking daily short-term stability of the instrument 9

11.1 General..... 9

11.2 Preparation and storage of sub-samples 9

11.3 Analysis of control samples 9

11.4 Monitoring the analytical procedure..... 9

11.5 Re-adjustment of calibration 11

12 Precision and accuracy..... 11

12.1 Repeatability..... 11

12.2 Intralaboratory reproducibility	12
12.3 Accuracy.....	12
13 Test report	12
Annex A (informative) Procedure for setting and evaluating linearity on a mass/volume basis	13
Annex B (informative) Control and adjustment of correction factors	15
Annex C (informative) Procedure for calibration of infrared milk analysers using modified milk samples	24
Bibliography.....	27

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 3.

Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

International Standard ISO 9622 was prepared by Technical Committee ISO/TC 34, *Agricultural food products*, Subcommittee SC 5, *Milk and milk products*, in collaboration with the International Dairy Federation (IDF) and AOAC International, and will also be published by these organizations.

Annexes A, B and C of this International Standard are for information only.

Whole milk — Determination of milkfat, protein and lactose content — Guidance on the operation of mid-infrared instruments

1 Scope

This International Standard describes the operating conditions for instruments used for the determination of fat, protein and lactose content of ex-farm milk, based upon the measurement of the absorption of mid-infrared radiation at wavelengths which are representative of each component analysed.

NOTE 1 In practice, these measurements are made using commercial automatic or semi-automatic instruments defined in clause 5 and called "infrared instruments" in this International Standard.

Any model of instrument which does not comply with the principle of analysis given in this International Standard or which incorporates modifications that may alter the principal characteristics of the apparatus (repeatability, accuracy, conditions of use), as well as the means of adjusting the calibration, will require a separate specific standard.

NOTE 2 Not all instruments allow determination of the lactose content. Moreover, some instruments allow direct measurement of the water content. The total solids content can be estimated by adding together the fat, protein and lactose contents, a constant being used to correct for the salt content.

The method described is applicable to the determination of the fat, protein and, as appropriate, lactose content of ex-farm milk. The method is also applicable to the analysis of milk of other species (goat, ewe, buffalo, etc.) and processed milk provided a specific calibration of the instrument (see clause 7) is made.

2 Normative references

The following normative documents contain provisions which, through reference in this text, constitute provisions of this International Standard. For dated references, subsequent amendments to, or revisions of, any of these publications do not apply. However, parties to agreement based on this International Standard are encouraged to investigate the possibility of applying the most recent editions of the normative documents indicated below. For undated references, the latest edition of the normative document referred to applies. Members of ISO and IEC maintain registers of currently valid International Standards.

ISO 1211, *Milk — Determination of fat content — Gravimetric method (Reference method)*.

ISO 5765-1, *Dried milk, dried ice-mixes and processed cheese — Determination of lactose content — Part 1: Enzymatic method utilizing the glucose moiety of the lactose*.

ISO 5765-2, *Dried milk, dried ice-mixes and processed cheese — Determination of lactose content — Part 2: Enzymatic method utilizing the galactose moiety of the lactose*.

ISO 8968-1, *Milk — Determination of nitrogen content — Part 1: Kjeldahl method*.

ISO 8968-2, *Milk — Determination of nitrogen content — Part 2: Block-digestion method (Macro method)*.

ISO 8968-4, *Milk — Determination of nitrogen content — Part 4: Determination of non-protein nitrogen content*.

ISO 8968-5, *Milk — Determination of nitrogen content — Part 5: Determination of protein-nitrogen content*.

3 Terms and definitions

For the purposes of this International Standard, the following terms and definitions apply.

3.1

infrared instrument

proprietary apparatus which, when used under the conditions defined in this International Standard, estimates the mass fraction of fat, protein and lactose in whole milk

3.2

fat, protein and lactose content

mass fraction of substances determined using the method specified in this International Standard

NOTE The fat, protein and lactose contents are expressed as mass fractions, in percent [formerly given as % (*m/m*)].

4 Principle

After homogenization of the milk sample, measurement with an infrared spectrometer of the quantity of radiation absorbed by:

- $\frac{3}{4}$ the carbonyl groups of the ester bonds of the glyceride at approximately 5,7 μm (traditionally referred to as the A Filter), and/or by the CH groups at approximately 3,5 μm (traditionally referred to as the B Filter), for determination of the fat content;
- $\frac{3}{4}$ the secondary amide groups of the peptide bonds at approximately 6,5 μm , for determination of the protein content;
- $\frac{3}{4}$ the hydroxyl groups of lactose at approximately 9,6 μm , for determination of the lactose content.

An estimate of the content of each component is made by reference to the amount of infrared light absorbed, either by water at the same wavelength or by the milk at a different wavelength at which there is only a slight absorption by the compound being measured.

NOTE For practical reasons, samples may be preserved with, for instance, a solution containing a mass fraction of 0,1 % potassium dichromate, 0,03 % sodium azide or 0,02 % to 0,06 % bronopol. It is necessary to check individual instrument response for all channels.

5 Principal characteristics of infrared instruments

The commercial instruments available may have one or two cells, with two wavebands or a single waveband per channel (component), and may use either a single- or a double-beam optical system, with either electronic ratioing or a servo system to estimate the transmitted radiation, and may produce the relevant wavelengths by diffraction grating, by interference optical filters, or by a Fourier-transformed interferogram. Instruments may differ also with respect to the number of specific wavelength ranges operating to predict the concentration of a given component.

NOTE In the case of interferometry, this International Standard is only applicable to the wavelength ranges mentioned in clause 4.

6 Factors affecting the accuracy of measurement

6.1 Instrument factors

6.1.1 Linearity

If an instrument is calibrated to express the results in mass/mass values, the solutions to adjust and evaluate linearity should be prepared on a mass/mass basis. If, on the other hand, the instrument is calibrated against volumetric reference methods or is calibrated to express results in mass/volume values, linearity should be set and evaluated on a mass/volume basis in order to give the optimal correlation to the reference method. Examples of how mass/volume evaluations can be carried out are given in annex A. This clause outlines only mass/mass based evaluations.

To check the linearity for each component, make up six solutions of known concentration as described in Table 1, using the following.

- Unhomogenized cream is recommended with a mass fraction of fat of 8 %, diluted with skimmed milk, to check the linearity at wavelengths of 5,7 μm (A Filter) and 3,5 μm (B Filter) for the determination of the fat content.
- UF skimmed milk retentate, diluted with ultrafiltrate is recommended to check the linearity at a wavelength of 6,5 μm for the determination of the protein content. Alternatively, whey protein concentrate, skim milk powder or evaporated skim milk diluted with distilled water may be used. The stock solution should contain a mass fraction of approx. 5,5 % protein.
- A solution of 60 g/l of lactose monohydrate, diluted with water, is recommended to check the linearity at a wavelength of 9,6 μm for the determination of the lactose content.

Table 1

Parts of stock solution (i.e. cream, UF retentate or 6 % lactose solution)	Parts of diluent (i.e. skim milk or water)	Relative concentration
100	0	1,0
80	20	0,8
60	40	0,6
40	60	0,4
20	80	0,2
0	100	0

The concentrations of the solutions should be in regular increments from zero to the desired upper limits of instrument readings.

Whenever possible, use the primary signal to check the linearity. Analyse each sample in triplicate and calculate the linear regression equation, as follows:

$$y = bx + a$$

and residuals r_i

$$r_i = y_i - (ax_i + b)$$

Plot residuals ϵ_i (y-axis) versus the concentration of the component in solution (x-axis) on a graph for each of the components. A visual inspection of the data points will usually yield sufficient information about the linearity of the signal.

If a more objective criterion for linearity is required, calculate the ratio of the residual range to the range of signal values:

$$\Delta\epsilon/\Delta s = (\epsilon_{\max} - \epsilon_{\min})/(s_{\max} - s_{\min})$$

where

$\epsilon_{\max}$ and $\epsilon_{\min}$ are the upper and lower residuals, respectively;

$s_{\max}$ and $s_{\min}$ are the upper and lower signal values, respectively.

Typical range for the ratio $\Delta\epsilon/\Delta s$ is 0,01 to 0,02.

Alternatively, a one-way analysis of variance can be carried out to confirm the non-linearity.

If the ratios between concentrations and instrument readings are not strictly linear over the whole range of measurements, adjust the linearity of the instrument response for the component, in accordance with the manufacturer's instructions.

NOTE Milk from animals other than cows may have higher concentrations of fat and protein. For such milks better performance may be obtained if the linearity is adjusted and checked specifically for the relevant concentration range.

6.1.2 Purging efficiency of the cell

After a single pumping sequence of a sample through the instrument cell, the residual volume of the previous samples shall not exceed 1 % of the total volume of the cell.

To check the effectiveness of rinsing, analyse 20 consecutive samples of water and homogenized whole milk, using the sequence: water, water, milk, milk, water, water, etc., and record for each sample of water and milk the readings at all wavelengths used. Calculate for each wavelength the purging efficiency, E , using the formula:

$$E = (\Sigma M_1 - \Sigma W_2)100/(\Sigma M_2 - \Sigma W_2)$$

where

M_1 is the first reading for milk;

M_2 is the second reading for milk;

W_2 is the second reading for water at the same wavelength.

The value of this ratio shall not be less than 99 %.

6.1.3 Homogenization

6.1.3.1 To check the efficiency of the homogenizer, make two consecutive analyses, first with an unhomogenized whole milk sample, and secondly with the same whole milk sample after it has been homogenized through the instrument's homogenizer. The difference between the two fat readings shall not exceed 0,05 % for a milk sample containing a mass fraction of 3,5 % milkfat. To calculate the appropriate pass/fail criteria for milkfat concentrations other than 3,5 %, multiply the actual fat content by 0,014 3 to obtain the new criteria.

NOTE This procedure is not applicable to some instruments.

CAUTION: The results of this test can be misleading, as an instrument in which the homogenizer does not work at all will give very little difference between the first and the second run. A safer but more laborious alternative method is described in 6.1.3.2 (see reference [2]).

6.1.3.2 Alternatively, obtain an unhomogenized as well as a homogenized portion of the same milk either by collecting raw and processed milk from the same tank at a dairy plant or by producing smaller volumes by means of a bench-top or pilot-plant homogenizer. Then measure both the unhomogenized and the same homogenized milk and compare the difference in results to the above-mentioned pass/fail criterion.

The assumption is that the homogenization efficiency of the external homogenizer is good. This can be verified by particle size analysis of the homogenized milk. An acceptable range for the diameter, d , of the fat globules in the reference homogenized milk is 0,75 μm to 0,85 μm , with a $d(0,9)$ of 1,4 μm to 1,5 μm [$d(0,9)$ means that 90 % of the milkfat has fat globules with a diameter of less than d]. At this particle size distribution, light scattering will be minimal at both the wavelengths corresponding to the determination of fat by filter A and by filter B.

6.1.4 Water vapour within the instrument

Variations in humidity of the air within the optical unit of the instrument result in variations in the optical zero and calibration. Replace the absorbent (silica gel) before it starts to change colour, preferably at regular intervals determined by trials. A once-a-week change, possibly before the week-end to allow time for the instrument to dry, is considered good practice.

6.2 Physicochemical and biological factors

6.2.1 Milk composition

The signal obtained at each wavelength is the result of the specific absorption by the component being determined and, to a greater or lesser extent, by the variations in concentration of the other major components, water included, and by salts.

The influence of the variations of the fat, protein and lactose content of milk is corrected by cross-correction or interaction factors that are specific to each wavelength and each type of instrument.

These interaction factors are calculated either by the manufacturer or by the user and are automatically introduced into the measurements.

The measurements to be carried out and the corresponding corrections that shall be applied are the following:

- ¾ no determination of the protein content at 6,5 μm without a simultaneous determination of the fat content and a correction of the protein reading for the fat content;
- ¾ no determination of the fat content at 3,5 μm without a simultaneous determination of the protein and lactose contents and a correction of the fat reading for the protein and lactose contents.

Although the adoption of other internal corrections is not considered compulsory, they are highly recommended because they significantly improve the accuracy of the measurements.

Every 3 months, check the correction factors of the instrument using, for instance, the methods described in annex B. The apparent interactions should be as close as possible to zero and should not exceed limits of $\pm 0,02$. Beyond these limits the cross-corrections shall be adjusted according to the manufacturer's recommendations.

The correction factors should be checked whenever any major part of the instrument, for instance the interference filters, is serviced or changed.

6.2.2 Fat

6.2.2.1 Fatty acid composition

The variations in the fatty acid composition of milk (mean molecular mass and degree of unsaturation) influence significantly the relationship between the results of the reference method and the infrared measurements at 5,7 μm and, to a lesser extent, at 3,5 μm .