

Status: Currently Official on 14-Feb-2025
 Official Date: Official as of 01-Feb-2023
 Document Type: USP Monographs
 DocId: GUID-0ABC20B3-883E-47DD-8456-2050BC79AC0E_4_en-US
 DOI: https://doi.org/10.31003/USPNF_M23490_04_01
 DOI Ref: g37yl

© 2025 USPC
 Do not distribute

Dexamethasone Sodium Phosphate Injection

DEFINITION

Dexamethasone Sodium Phosphate Injection is a sterile solution of Dexamethasone Sodium Phosphate in Water for Injection. It contains NLT 90.0% and NMT 115.0% of the labeled amount of dexamethasone phosphate ($C_{22}H_{30}FO_8P$), present as the disodium salt.

IDENTIFICATION

• A. THIN-LAYER CHROMATOGRAPHY

Solution A: Mix 3.1 g of boric acid and 500 mL of water in a 1-L volumetric flask. Add 21 mL of 1 N sodium hydroxide and 10 mL of 0.1 M magnesium chloride, and dilute with water to volume.

Solution B: 1 mg/mL of alkaline phosphatase enzyme in *Solution A*

Standard solution: 300 µg/mL of [USP Dexamethasone RS](#) in methylene chloride

Sample solution: Pipet a volume of Injection, equivalent to 10 mg of dexamethasone phosphate, into a 100-mL volumetric flask, and dilute with water to volume. Pipet 5 mL of this solution into a 125-mL separator, and wash with two 10-mL portions of water-washed methylene chloride, discarding the washings. Transfer the solution to a glass-stoppered, 50-mL tube, and add 5 mL of *Solution B*. Allow to stand at 37° for 45 min, and extract with 25 mL of methylene chloride. Evaporate 15 mL of the methylene chloride extract on a steam bath to dryness, and dissolve the residue in 1 mL of methylene chloride.

Chromatographic system

(See [Chromatography \(621\)](#), [Thin-Layer Chromatography](#).)

Adsorbent: 0.25-mm layer of chromatographic silica gel mixture (20- × 20-cm plate)

Application volume: 5 µL

Developing solvent system: Chloroform, acetone, and water (50:50:1)

Spray reagent: Dilute sulfuric acid (1 in 2)

Analysis

Samples: *Standard solution* and *Sample solution*

Allow the spots to dry, and develop the chromatogram using the *Developing solvent system* in a tank completely lined with filter paper, until the solvent front has moved about three-fourths of the length of the plate. Remove the plate from the developing tank, mark the solvent front, and allow the spots to dry. Spray the plate with *Spray reagent*, and heat at 105° until brown or black spots appear.

Acceptance criteria: The R_f value of the principal spot of the *Sample solution* corresponds to that of the *Standard solution*.

• **B.** The retention time of the major peak of the *Sample solution* corresponds to that of the *Standard solution*, as obtained in the Assay.

ASSAY

Change to read:

• PROCEDURE

Mobile phase: 0.01 M monobasic potassium phosphate in a mixture of methanol and water (1:1)

Standard solution: 0.09 mg/mL of freshly prepared [USP Dexamethasone Sodium Phosphate RS](#) in *Mobile phase*

Sample solution: Nominally 0.08 mg/mL of dexamethasone phosphate from Injection, in *Mobile phase*

Chromatographic system

(See [Chromatography \(621\)](#), [System Suitability](#).)

Mode: LC

Detector: UV 254 nm

Column: 4-mm × 30-cm; packing L1

Flow rate: 1.6 mL/min

Injection volume: 20 µL

System suitability

Sample: *Standard solution*

[NOTE—The retention time for dexamethasone phosphate is about 5 min.]

Suitability requirements

Relative standard deviation: NMT 1.5%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the percentage of the labeled amount of dexamethasone phosphate ($C_{22}H_{30}FO_8P$) in the portion of Injection taken:

Result = $(r_U/r_S) \times (C_S/C_U) \times (M_{r1}/M_{r2}) \times 100$ (ERR 1-Feb-2023)

- r_U = peak response from the *Sample solution*
- r_S = peak response from the *Standard solution*
- C_S = concentration of [USP Dexamethasone Sodium Phosphate RS](#) in the *Standard solution* (mg/mL)
- C_U = nominal concentration of dexamethasone phosphate in the *Sample solution* (mg/mL)
- M_{r1} = molecular weight of dexamethasone phosphate, 472.44
- M_{r2} = molecular weight of dexamethasone sodium phosphate, 516.40

Acceptance criteria: 90.0%–115.0%

SPECIFIC TESTS

- pH** (791): 7.0–8.5
- BACTERIAL ENDOTOXINS TEST** (85): NMT 31.3 USP Endotoxin Units/mg of dexamethasone phosphate
- OTHER REQUIREMENTS:** Meets the requirements in [Injections and Implanted Drug Products \(1\)](#).

ADDITIONAL REQUIREMENTS

- PACKAGING AND STORAGE:** Preserve in single-dose or multiple-dose containers, preferably of Type I glass, protected from light. Store at controlled room temperature.
- USP REFERENCE STANDARDS** (11):
 - [USP Dexamethasone RS](#)
 - [USP Dexamethasone Sodium Phosphate RS](#)

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
DEXAMETHASONE SODIUM PHOSPHATE INJECTION	Documentary Standards Support	SM52020 Small Molecules 5

Chromatographic Database Information: [Chromatographic Database](#)

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. PF 43(3)

Current DocID: GUID-0ABC20B3-883E-47DD-8456-2050BC79AC0E_4_en-US

DOI: https://doi.org/10.31003/USPNF_M23490_04_01

DOI ref: [g37yl](#)