

INTERNATIONAL STANDARD

**ISO
15163**

**IDF
110**

First edition
2012-05-15

Milk and milk products — Calf rennet and adult bovine rennet — Determination by chromatography of chymosin and bovine pepsin contents

*Lait et produits laitiers — Présure de veau et coagulant issu de bovin
adulte — Détermination des teneurs en chymosine et en pepsine
bovine par chromatographie*



Reference numbers
ISO 15163:2012(E)
IDF 110:2012(E)

AVA **Documents**
© ISO and IDF 2012
Document Sharing Hub



COPYRIGHT PROTECTED DOCUMENT

© ISO and IDF 2012

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from either ISO or IDF at the respective address below.

ISO copyright office
Case postale 56 • CH-1211 Geneva 20
Tel. + 41 22 749 01 11
Fax + 41 22 749 09 47
E-mail copyright@iso.org
Web www.iso.org

Published in Switzerland

International Dairy Federation
Silver Building • Boulevard Auguste Reyers 70/B • B-1030 Brussels
Tel. + 32 2 733 98 88
Fax + 32 2 733 04 13
E-mail info@fil-idf.org
Web www.fil-idf.org

Contents

Page

Foreword	iv
Foreword	v
Introduction	vi
1 Scope	1
2 Normative references	1
3 Principle	1
4 Reagents	2
5 Apparatus	3
6 Sampling	4
7 Procedure	4
7.1 Check	4
7.2 Preparation of a fresh column with Fractogel	4
7.3 Regeneration and equilibration of the Fractogel resin in the column	4
7.4 Storage of the column with Fractogel	5
7.5 Preparation of test sample	5
7.6 Analysis of the desalted rennet	6
7.7 Determination of the clotting times	7
8 Calculation and expression of results	8
8.1 Calculation of the activity of chymosin and pepsin, expressed as a percentage	8
8.2 Calculation of active chymosin and active bovine pepsin, in milligrams per litre	9
8.3 Expression of results	10
9 Precision	10
9.1 Interlaboratory test	10
9.2 Repeatability	10
9.3 Reproducibility	11
10 Test report	11
Annex A (informative) Qualitative determination by double immunodiffusion of milk-coagulating enzymes in commercial coagulants	12
Annex B (informative) Interlaboratory test	17
Bibliography	19

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO 15163|IDF 110 was prepared by Technical Committee ISO/TC 34, *Food products*, Subcommittee SC 5, *Milk and milk products* and the International Dairy Federation (IDF). It is being published jointly by IDF and ISO.

This first edition of ISO 15163|IDF 110 cancels and replaces IDF 110B:1997, which has been technically revised.

Foreword

IDF (the International Dairy Federation) is a non-profit organization representing the dairy sector worldwide. IDF membership comprises National Committees in every member country as well as regional dairy associations having signed a formal agreement on cooperation with IDF. All members of IDF have the right to be represented on the IDF Standing Committees carrying out the technical work. IDF collaborates with ISO in the development of standard methods of analysis and sampling for milk and milk products.

The main task of Standing Committees is to prepare International Standards. Draft International Standards adopted by the Standing Committees are circulated to the National Committees for endorsement prior to publication as an International Standard. Publication as an International Standard requires approval by at least 50 % of IDF National Committees casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. IDF shall not be held responsible for identifying any or all such patent rights.

ISO 15163|IDF 110 was prepared by the International Dairy Federation (IDF) and Technical Committee ISO/TC 34, *Food products*, Subcommittee SC 5, *Milk and milk products*. It is being published jointly by IDF and ISO.

All work was carried out by an ISO-IDF Project Group on *Chymosin and bovine pepsin determination*, of the Standing Committee on *Analytical methods for processing aids and indicators*, under the aegis of its project leader, Prof. A. Andrén (SE).

This first edition of ISO 15163|IDF 110 cancels and replaces IDF 110B:1997, which has been technically revised.

Introduction

Calf rennet and adult bovine rennet preparations contain both chymosin and bovine pepsin in various amounts as main clotting enzymes. The proportion of chymosin decreases relative to pepsin in the abomasum (the fourth “true” stomach) with age and at weaning of the calf.

The ratio of abomasa from young cattle to that of old cattle in the raw material for rennet production thus highly influences the composition of chymosin and pepsin in the final rennet. The higher the abomasa from young milk-fed calves, the higher the proportion of chymosin and vice versa^{[5][6]}.

Both chymosin and pepsin have special characteristics relevant to milk-clotting activity and suitability for cheese making. The milk-clotting activity of pepsin is, for example, much more pH-dependent than chymosin and pepsin also has a more general proteolytic activity than chymosin.

Therefore, it is very important to analyse the content of chymosin and pepsin in addition to the strength (total milk-clotting activity) of the rennet^{[6][7]}.

Milk and milk products — Calf rennet and adult bovine rennet — Determination by chromatography of chymosin and bovine pepsin contents

1 Scope

This International Standard specifies a reference method for the determination of the amounts of chymosin and bovine pepsin present in a test sample of calf rennet and adult bovine rennet. In addition, it can be used for mixtures of calf/bovine rennet with fermentation-produced bovine chymosin (FPC).

2 Normative references

The following referenced document is indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 11815|IDF 157:2007, *Milk — Determination of total milk-clotting activity of bovine rennets*

3 Principle

As a first step, the rennet sample is desalted and the enzymes chymosin and bovine pepsin separated on an anion exchange column^{[8][9]}. In a second step, the milk-clotting activity of each of the separated two enzymes is determined by ISO 11815|IDF 157 (reconstituted milk with pH 6,5). The enzymatic composition of the rennet sample is expressed in percentage chymosin activity and percentage pepsin activity of the sum of the activities in International Milk-Clotting Units (IMCU) of both components, or the results are expressed in milligrams per litre of active chymosin and milligrams per litre of active pepsin.

The total milk-clotting activity of the first batch of calf rennet reference standard powder and the first batch of adult bovine rennet reference standard powder has once and for all been set at 1 000 IMCU/g. Future preparations of reference standards shall be set relative to the previous reference standards (see ISO 11815|IDF 157).

This International Standard specifies both a manual set-up of the anion exchange chromatography and an alternative automated set-up.

This is a reference method and changes may therefore only be made if confirmed to give the same result, and repeatability and reproducibility at least as high as the original standard method. Any change to what is stated in this International Standard method shall also be mentioned in the test report (see Clause 10).

4 Reagents

Use only reagents of recognized analytical grade, unless otherwise specified, and distilled or demineralized water or water of equivalent purity.

4.1 Resin, Fractogel® EMD DEAE (M) (Merck cat. no. 1.16883)¹⁾ or Mono Q® 1 ml prepacked column (HR 5/5 or 5/50 GL from GE Healthcare)²⁾ or equivalent resin.

NOTE 1 Fractogel® EMD DEAE (M) is a suitable resin for manual chromatography and Mono Q® is suitable for the automated chromatography.

NOTE 2 If the Fractogel® or Mono Q resins are substituted by another resin, it will most likely be necessary to change the buffers in 4.12, resulting in a need to re-evaluate the method.

4.2 Piperazine hexahydrate ($C_4H_{10}N_2 \cdot 6H_2O$).

4.3 Sodium chloride (NaCl).

4.4 Thymol, optional preservative.

4.5 Sodium hydroxide (NaOH).

4.6 Hydrochloric acid solution, $c(HCl) = 1 \text{ mol/l}$.

4.7 Ethanol (C_2H_5OH), with a volume fraction of at least 96 %.

4.8 Ethanol (C_2H_5OH), with an approximate volume fraction of at least 20 %.

Add 105 ml ethanol 96 % (4.7) to 400 ml water and mix. If a sterile filtration is desired, filter the water before mixing it with the ethanol.

4.9 Urea, $c(N_2H_4CO) = 8 \text{ mol/l}$.

Dissolve 48 g urea in water and fill to a total volume of 100 ml.

4.10 Dialysis tubing, of diameter approximately 1 cm (Union Carbide)³⁾ or equivalent (optional).

NOTE The quality of the dialysis tubing is not critical.

4.11 Desalting columns, Bio-Rad – Econopac 10DG (cat. no. 732-2010)⁴⁾ or equivalent (optional).

Use either the dialysis tubing (4.10) or the desalting columns for desalting the rennet.

4.12 Buffer solutions

4.12.1 Buffer solution I, piperazine [$(CH_2)_4(NH)_2$], $c[(CH_2)_4(NH)_2] = 0,025 \text{ mol/l}$.

1) Fractogel® EMD DEAE (M) is an example of a suitable product available commercially. This information is given for the convenience of the users of this document and does not constitute an endorsement of the product by ISO or IDF.

2) Mono Q® 1 ml prepacked column is an example of a suitable product available commercially. This information is given for the convenience of the users of this document and does not constitute an endorsement of the product by ISO or IDF.

3) Union Carbide is an example of a suitable product available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement of the product by ISO and IDF.

4) Bio-RAD - Econopac 10DG is an example of a suitable product available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement of the product by ISO or IDF.