



JAPANESE
INDUSTRIAL
STANDARD

Translated and Published by
Japanese Standards Association

JIS A 1966 : 2005

**Indoor air—Sampling and analysis
of volatile organic compounds by
sorbent tube/thermal desorption /
capillary gas chromatography—
Pumped sampling**

ICS 13.040.20

Reference number : JIS A 1966 : 2005 (E)

PROTECTED BY COPYRIGHT

 **Documents**
Document Sharing Hub

24 5

Foreword

This translation has been made based on the original Japanese Industrial Standard revised by the Minister of Land, Infrastructure and Transportation through deliberations at the Japanese Industrial Standards Committee.

This Standard has been made based on **ISO 16017-1 : 2000** *Indoor, ambient and workplace air—Sampling and analysis of volatile organic compounds by sorbent tube/thermal desorption/thermal desorption/capillary gas chromatography—Part 1: Pumped sampling* for the purpose of making it easier to compare this Standard with International Standard; to prepare Japanese Industrial Standard conforming with International Standard; and to propose a draft of an International Standard which is based on Japanese Industrial Standard.

Attention is drawn to the possibility that some parts of this Standard may conflict with a patent right, application for a patent after opening to the public, utility model right or application for registration of utility model after opening to the public which have technical properties. The relevant Minister and the Japanese Industrial Standards Committee are not responsible for identifying the patent right, application for registration of utility model after opening to the public which have the said technical properties.

Date of Establishment: 2005-11-01

Date of Public Notice in Official Gazette: 2005-11-01

Investigated by: Japanese Industrial Standards Committee
Standards Board
Technical Committee on Architecture

JIS A 1966 : 2005, First English edition published in 2006-09

Translated and published by: Japanese Standards Association
4-1-24, Akasaka, Minato-ku, Tokyo, 107-8440 JAPAN

In the event of any doubts arising as to the contents,
the original JIS is to be the final authority.

© JSA 2006

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from the publisher.

Printed in Japan

NH/HN

PROTECTED BY COPYRIGHT

Contents

	Page
Introduction	1
1 Scope	1
2 Normative references	2
3 Terms and definitions	2
4 Principle	3
5 Reagents and materials	3
5.1 Volatile organic compounds (VOC)	3
5.2 Dilution solvent	3
5.3 Sorbents	3
5.4 Calibration standards	4
5.5 Standard atmosphere	4
5.6 Standard atmosphere loaded sorbent tubes	4
5.7 Preparation of solutions for liquid spiking	4
5.8 Standard liquid loaded sorbent tubes	5
6 Apparatus and instruments	6
6.1 Sorbent tubes	6
6.2 Sorbent tube end caps	6
6.3 Sorbent tube unions	6
6.4 Syringes	6
6.5 Sampling pump	6
6.6 Plastic or rubber tubing	6
6.7 Soap-bubble meter or other suitable device for calibrating pump	6
6.8 Gas chromatograph	7
6.9 Thermal desorption apparatus	7
6.10 Injection facility for preparing standard tube	7
7 Tube conditioning	7
8 Calibration of pump	7
9 Sampling	7
10 Procedure	9
10.1 Safety precautions	9
10.2 Desorption and analysis	9
10.3 Calibration	11
10.4 Determination of sample concentration	12
10.5 Determination of desorption efficiency	12

11	Calculations	12
11.1	Mass concentration of analyte	12
11.2	Volume concentration of analyte	13
12	Interferences	13
13	Performance characteristics	14
14	Test report	14
15	Quality control	14
Annex A (normative)	Determination of breakthrough volumes from standard atmosphere	26
Annex B (normative)	Presumption of retention volume by extrapolation	27
Annex C (informative)	Description of sorbent types	28
Annex D (informative)	Guidance on sorbent selection	29
Annex E (informative)	Guidance on sorbent use	30
Annex F (informative)	Summary on overall uncertainty, precision, bias and storage	32
Bibliography.....		35
Annex 1 (informative)	Comparison table between JIS and corresponding International Standard	37

Indoor air—Sampling and analysis of volatile organic compounds by sorbent tube/thermal desorption/capillary gas chromatography—Pumped sampling

Introduction This Japanese Industrial Standard has been prepared based on the first edition of **ISO 16017-1** *Indoor, ambient and workplace air—Sampling and analysis of volatile organic compounds by sorbent tube/thermal desorption/capillary gas chromatograph—Part 1: Pumped sampling* published in 2000 with some modifications of the technical contents.

The portions given dotted underlines are the matters in which the contents of the original International Standard have been modified. A list of modifications with the explanations is given in Annex 1 (informative).

1 Scope This Standard gives general guidance for the sampling and analysis of volatile organic compounds (VOCs) in air and shall be applied to indoor air.

This Standard is appropriate for a wide range of VOCs, including hydrocarbons, halogenated hydrocarbons, esters, glycol ethers, ketones and alcohols. Sorbent tubes having different applicability ranges are recommended for the sampling of these VOCs. Very low boiling compounds will only be partially retained by the sorbents and can only be estimated qualitatively. Semi-volatile compounds will be fully retained by the sorbents, but may only be partially recovered. Compounds for which this Standard has been tested are given in tables 1 to 13. This Standard may be applicable to compounds not listed, but in these cases it is advisable to use a back-up tube containing the same or a stronger sorbent.

This Standard is applicable to the measurement of airborne vapours of VOCs in a concentration range of approximately $0.5 \mu\text{g}/\text{m}^3$ to $10^5 \mu\text{g}/\text{m}^3$ individual compound.

The upper limit of the useful range is set by the sorptive capacity of the sorbent used and by the linear dynamic range of the gas chromatograph column and detector or by the sample-splitting capability of the analytical instrumentation used. The sorptive capacity is measured as a breakthrough volume of air, which determines the maximum air volume that shall not be exceeded when sampling.

The lower limit of the useful range depends on the noise level of the detector and interfering artefacts from the analysis system or the sorbent tubes. Artefacts are typically 1 ng for well-conditioned Tenax GR and carbonaceous sorbents such as Carboxpack/Carbotrap type materials, carbonized molecular sieves and molecular sieves such as Spherosorb, or pure charcoal; at a few ng levels for Tenax TA and at 5 ng to 50 ng levels for other porous polymers such as Chromosorbs and Porapak.

The procedure specified in this Standard is applicable to low flowrate personal sampling pumps and gives a time-weighted average result. It is not applicable to the measurement of instantaneous or short-term fluctuations in concentration.

- Remarks 1 The sorbents listed in Annex C and elsewhere in this Standard are those known to perform as specified under this Standard. Each sorbent or product that is identified by a trademarked name is unique and has a sole manufacturer, however, they are widely available from many different suppliers. This information is given for the convenience of users of this Standard and does not constitute an endorsement of the product named. Equivalent products may be used if they can be shown to lead to the same results.
- 2 The International Standard corresponding to this Standard is as follows.

In addition, symbols which denote the degree of correspondence in the contents between the relevant International Standard and **JIS** are IDT (identical), MOD (modified), and NEQ (not equivalent) according to **ISO/IEC Guide 21**.

ISO 16017-1:2000 *Indoor, ambient and workplace air—Sampling and analysis of volatile organic compounds by sorbent tube /thermal desorption /capillary gas chromatography—Part 1: Pumped sampling (MOD)*

2 Normative references The following standards contain provisions which, through reference in this text, constitute provisions of this Standard. Only the edition of the indicated year constitutes the provision of this Standard but the revision and amendment made thereafter are not apply.

- ISO 6141:2000** *Gas analysis—Requirements for certificates for calibration gases and gas mixtures*
- ISO 6145-1:1986** *Gas analysis—Preparation of calibration gas mixtures—Dynamic volumetric methods—Part 1: Methods of calibration*
- ISO 6145-3:1986** *Gas analysis—Preparation of calibration gas mixtures—Dynamic volumetric methods—Part 3: Periodic injections into a flowing gas stream*
- ISO 6145-4:1986** *Gas analysis—Preparation of calibration gas mixtures—Dynamic volumetric methods—Part 4: Continuous injection method*
- ISO 6145-5:1986** *Gas analysis—Preparation of calibration gas mixtures—Dynamic volumetric methods—Part 5: Capillary calibration devices*
- ISO 6145-6:1986** *Gas analysis—Preparation of calibration gas mixtures—Dynamic volumetric methods—Part 6: Critical orifices*
- ISO 6349:1979** *Gas analysis—Preparation of calibration gas mixtures—Permeation method*
- EN 1076:1997** *Workplace atmospheres—Pumped sorbent tubes for the determination of gases and vapours—Requirements and test methods*

3 Terms and definitions For the purposes of this Standard, the following definitions apply.

3.1 breakthrough volume air volume that passed through a sorbent tube before the concentration of eluting vapour reaches 5 % of the applied sample concentration

Remarks 1 The breakthrough volume varies with the vapour and the sorbent type.

2 See bibliography [4].

3.2 retention volume carrier gas volume at peak maximum of a small aliquot of an organic vapour eluted from a sorbent tube by chromatographic carrier gas

3.3 Safe Sampling Volume (SSV) the maximum sampling volume of a sorbent tube to be used, or the sampled air volume equivalent to 70 % of breakthrough volume or 50 % of retention volume

4 Principle A measured volume of sample air is drawn through one (or more) sorbent tubes in series; an appropriate sorbent (or sorbents) is selected for the compound or mixture to be sampled. Provided suitable sorbents are chosen, volatile organic compounds (VOC) are retained by the sorbent tube and thus are removed from the air stream. The collected vapour on each tube is desorbed by heat and is transferred under inert carrier gas into a gas chromatograph equipped with a capillary column and a flame ionization detector or other suitable detector, where it is analysed. Calibration is achieved by means of liquid or vapour spiking onto a sorbent tube.

5 Reagents and materials During the analysis, use only reagents of recognized analytical reagent grade. Calibration blend solutions should be prepared weekly, or more frequently if evidence is noted of deterioration, such as condensation reactions between alcohols and carboxylic acid.

5.1 Volatile organic compounds (VOC) VOC are required for calibration reagents, using either standard liquid spiking (5.7 and 5.8) or standard atmosphere spiking (5.4 to 5.6) onto sorbent tubes.

5.2 Dilution solvent Dilution solvent shall be used for preparing calibration blend solution for standard liquid spiking (5.7). This should be of chromatographic quality. It shall be free from compounds co-eluting with the analyte (5.1).

Remarks : Methanol is usually used. Alternative dilution solvents can be used, particularly if there is no possibility of reaction or chromatographic co-elution.

5.3 Sorbents Sorbents of particle size 0.18 to 0.25 mm (60 to 80 mesh) are recommended. Each sorbent shall be preconditioned under a flow of inert gas by heating it overnight at a temperature at least 25 °C below the published maximum for that sorbent before packing the tubes. They shall be kept in a clean atmosphere during cooling to room temperature, loaded into the sorbent tube, then stored. Wherever possible, analytical desorption temperatures should be kept below those of aging condition.

Tubes prepacked by the manufacturer are also available and as such only require conditioning.

Remarks 1 Sorbent particle sizes larger than 0.18 mm to 0.25 mm may be used but the breakthrough characteristics given in tables 1 to 6 may be affected. Smaller sorbent particle size ranges are not recommended because of back-pressure problems.

- 2 A description of sorbents is given in Annex C and a guide for sorbent selection is given in Annex D. Equivalent sorbents may be used. A guide to sorbent conditioning and analytical desorption parameters is given in Annex E.

5.4 Calibration standards Calibration standards are preferably prepared by loading required amounts of the analyte on sorbent tubes from standard atmospheres (see **5.5** and **5.6**), as this procedure most closely resembles the practical sampling situation.

If this way of preparation is not practicable, standards may be prepared by a liquid spiking procedure (see **5.7** and **5.8**), provided that the accuracy is confirmed by either:

- a) the method established by using procedures giving spiking levels fully traceable to primary standards of mass and/or volume
- b) the method confirmed by reference materials
- c) the method confirmed by standards produced using standard atmospheres
- d) the method confirmed by results of reference measurement procedures

Remarks : The loading ranges given in **5.6**, **5.7** and **5.8** are not mandatory. For specific applications in which larger volumes are used to measure lower concentrations, other loading ranges may be more appropriate.

5.5 Standard atmosphere Prepare standard atmosphere of known concentrations of the analyte by a recognized procedure. Methods described in **ISO 6141**, **ISO 6145** and **ISO 6349** are suitable. If the procedure is not applied under conditions that allow the establishment of full traceability of the generated concentrations to primary standards of mass and/or volume, or if the chemical inertness of the generation system cannot be guaranteed, the concentrations shall be confirmed using an independent procedure.

5.6 Standard atmosphere loaded sorbent tubes Prepare standard atmosphere loaded sorbent tubes by passing an accurately known volume of the calibration atmosphere through the sorbent tube, e.g. by means of a pump. The volume of atmosphere sampled shall not exceed the breakthrough volume of the sorbent. After loading, disconnect and seal the tube. Prepare fresh standards with each batch of samples.

For example, prepare standard atmosphere equivalent to 10 mg/m³ and 100 µg/m³ of the analyte. For workplace air, load sorbent tubes with 100 ml, 200 ml, 400 ml, 1 L, 2 L, or 4 L of the 10 mg/m³ atmosphere. For ambient or indoor air, load sorbent tubes with 100 ml, 200 ml, 400 ml, 1 L, 2 L, 4 L or 10 L of the 100 µg/m³ atmosphere.

5.7 Preparation of solutions for liquid spiking

5.7.1 Solution containing approximately 10 mg/ml of each liquid component Accurately weigh approximately 1 g of a substance or substances of interest into a 100 ml volumetric flask, starting with the least volatile substance. Make up to 100 ml with dilution solvent (**5.2**), stopper and shake to mix.

5.7.2 Solution containing approximately 1 mg/ml of each liquid component Introduce 50 ml of dilution solvent into a 100 ml volumetric flask. Add 10 ml of solution **5.7.1**. Make up to 100 ml with dilution solvent, stopper and shake to mix.

5.7.3 Solution containing approximately 100 µg/ml of each liquid component

Accurately weigh approximately 10 mg of a substance or substances of interest into a 100 ml volumetric flask, starting with the least volatile substance. Make up to 100 ml with dilution solvent (5.2), stopper and shake to mix.

5.7.4 Solution containing approximately 10 µg/ml of each liquid component

Introduce 50 ml of dilution solvent into a 100 ml volumetric flask. Add 10 ml of solution described in 5.7.3. Make up to 100 ml with dilution solvent, stopper and shake to mix.

5.7.5 Solution containing approximately 1 mg/ml of gas components

For gases, e.g. ethylene oxide, a high-level calibration solution may be prepared as follows. Obtain pure gas at atmospheric pressure by filling a small plastic gas bag from a gas cylinder. Fill a 1 mL gas-tight syringe with 1 ml of the pure gas and close the valve of the syringe. Using a 2 ml septum vial, add 2 ml dilution solvent and close with the septum cap. Insert the tip of the syringe needle through the septum cap into the dilution solvent. Open the valve and withdraw the plunger slightly to allow the dilution solvent to enter the syringe. The action of the gas dissolving in the dilution solvent creates a vacuum, and the syringe fills with solvent. Return the solution to the vial. Flush the syringe twice with the solution and return the washings to the vial. Calculate the mass of gas added using the gas laws, i.e. 1 mole of gas at STP (standard temperature and pressure: 273.15 K and 1 013.25 hPa) occupies 22.4 L. Correct for any non-ideality of the particular pure gas compound.

5.7.6 Solution containing approximately 10 µg/ml of gas components

For gases, e.g. ethylene oxide, a low-level calibration solution may be prepared as follows. Obtain pure gas at atmospheric pressure by filling a small plastic gas bag from a gas cylinder. Fill a 10 µl gas-tight syringe with 1 µl of the pure gas and close the valve of the syringe. Using a 2 ml septum vial, add 2 ml dilution solvent and close with the septum cap. Insert the tip of the syringe needle through the septum cap into the dilution solvent. Open the valve and withdraw the plunger slightly to allow the dilution solvent to enter the syringe. The action of the gas dissolving in the dilution solvent creates a vacuum, and the syringe fills with solvent. Return the solution to the vial. Flush the syringe twice with the solution and return the washings to the vial. Calculate the mass of gas added using the gas laws, i.e. 1 mole of gas at STP occupies 22.4 L. Correct for any non-ideality of the particular pure gas compound.

5.8 Standard liquid loaded sorbent tubes

Standard liquid loaded sorbent tubes are prepared by injecting aliquots of standard solutions onto clean sorbent tubes as follows. A sorbent tube is fitted into the injection unit (6.10) through which inert purge gas and a 1 µl to 4 µl aliquot of an appropriate standard solution, injected through the septum, are passed. After an appropriate time, disconnect and seal the tube. Prepare fresh standard tube with each batch of samples. For workplace air, load sorbent tubes with 1 µl to 5 µl of solutions 5.7.1, 5.7.2 or 5.7.5. For ambient and indoor air, load sorbent tubes with 1 µl to 5 µl of solutions 5.7.3, 5.7.4 or 5.7.6.

Remarks : In the case that the dilution solvent is methanol, a purge gas flowrate of 100 ml/min and a 5 min purge time have been found to be appropriate to eliminate most of the solution solvent from the tube. If other dilution solvents are used, the purge conditions should be determined experimentally.

6 Apparatus and instruments Use the following laboratory instruments and apparatus.

6.1 Sorbent tubes Typically, but not exclusively, sorbent tubes are constructed of stainless steel tubing 6.3 mm outside diameter, 5 mm inside diameter and 90 mm long. Tubes of other dimensions may be used but the safe sampling volumes (SSV) given in tables 1 to 6 are based on these tube dimensions. For labile analytes, such as unstable compounds, glass-lined or glass tubes (typically 4 mm inside diameter) should be used. One end of the tube is marked, for example by a scored ring about 10 mm from the sampling inlet end. The tubes are packed with preconditioned sorbents (5.3) so that the sorbent bed will be within the desorber heated zone and a gap of at least 14 mm is retained at each end to minimize errors due to diffusive ingress at very low pump flowrates.

Tubes contain 200 to 1 000 mg sorbent, depending on sorbent density (typically about 250 mg porous polymer or 500 mg carbon molecular sieve or graphitized carbon). The sorbents are retained by stainless steel gauzes and/or unsilanized glass wool plugs. If more than one sorbent is used in a single tube, the sorbents should be arranged in order of increasing sorbent strength and separated by unsilanized glass wool, with the weakest sorbent nearest to the marked sampling inlet end.

Do not pack sorbents with widely different ($>50\text{ }^{\circ}\text{C}$) maximum desorption temperatures into a single tube, or it will be impossible to condition or desorb the more stable sorbent(s) sufficiently thoroughly without causing degradation of the least stable sorbent(s).

6.2 Sorbent tube end caps The tubes shall be sealed according to the requirements of **EN 1076**, or equivalent, e.g. with metal screw-cap fittings with polytetrafluoroethylene (PTFE) seals.

6.3 Sorbent tube unions Two sorbent tubes may be connected in series during sampling with metal screw-cap couplings with PTFE seals.

6.4 Syringes Use syringes including a precision 10 μl liquid syringe readable to 0.1 μl , a precision 10 μl gas-tight syringe readable to 0.1 μl and a precision 1 ml gas-tight syringe readable to 0.01 ml.

6.5 Sampling pump The pump should fulfil the requirements of **EN 1232** [10] or equivalent.

6.6 Plastic or rubber tubing Plastic or rubber tubing is about 90 cm long, of appropriate diameter to ensure a leak-proof fit to both pump and tube or tube holder, if used. Clips should be provided to hold the tube and connecting tubing.

Tubes shall not be used with plastic or rubber tubing upstream of the sorbent. The use of such tubing may introduce contaminants or sorbed sampled VOCs.

6.7 Soap-bubble meter or other suitable device for calibrating pump The flow meter shall be traceably calibrated to a primary flow standard.

Remarks : The use of an uncalibrated integral flow meter for the calibration of pump flowrates may result in systematic errors of several tens of percent.

6.8 Gas chromatograph A gas chromatograph shall be fitted with a flame ionization, photoionization detector, mass spectrometric or other suitable detector, capable of detecting an injection of 0.5 ng toluene with an S/N ratio of at least 5 to 1.

The gas chromatograph shall have a capillary column capable of separating the analyte from other substances.

6.9 Thermal desorption apparatus Thermal desorption apparatus is for the two-stage thermal desorption of the sorbent tubes and for transfer of the desorbed vapours via an inert gas flow into a gas chromatograph. A typical apparatus contains a mechanism for holding the tubes to be desorbed whilst they are heated and purged simultaneously with inert carrier gas. The desorption temperature and time is adjustable, as is the carrier gas flowrate. The apparatus should also incorporate additional features, such as automatic sample tube loading, leak testing, and a cold trap in the transfer line to concentrate the desorbed sample. The desorbed sample, contained in the purge gas, is routed to the gas chromatograph and capillary column via a heated transfer line.

6.10 Injection facility for preparing standard tube A gas chromatographic injection port may be used for preparing standard tube. This can be used *in situ*, or it can be mounted separately. The carrier gas line to the injector should be retained. The back of the injection facility should be adapted if necessary to fit the tube. This can be done conveniently by means of a compression coupling with an O-ring seal.

7 Tube conditioning Prior to use, tubes should be conditioned by desorbing them at a temperature at or just above the analytical desorption temperature (see Annex E). For example, they should be reconditioned by the carrier gas like of 100 ml/min flowrate for 10 min. The carrier gas flow should be in a direction to the sampling end to prevent the recontamination of sorbents.

Tubes should then be analysed, using routine analytical parameters, to ensure that the thermal desorption blank is sufficiently small. If the blank is unacceptable, tubes should be reconditioned by repeating this procedure.

Once a sample has been analysed, the tube may be reused to collect a further sample immediately. However, it is advisable to check the thermal desorption blank if the tubes are left for an extended period before reuse, or if sampling for a different analyte is envisaged. Tubes should be sealed with metal screwcaps with combined PTFE ferrule fittings and stored in an airtight container when not used for sampling or being conditioned.

Remarks : The sorbent tube blank level is acceptable if interfering peaks are no greater than 10 % of the typical areas of the analytes of interest.

8 Calibration of pump Calibrate the pump with a representative sorbent tube assembly in-line, using an appropriate external calibrated meter. One end of the calibrated flow meter should be at atmospheric pressure to ensure proper operation.

9 Sampling Select a sorbent tube (or tube combination) appropriate for the compound or mixture to be sampled. Guidance on suitable sorbents is given in Annex D.

If more than one tube is to be used, prepare a tube assembly by joining the tubes with a union (**6.3**). Attach the pump to the sorbent tube or tube assembly with plastic or rubber tubing, so that the tube containing the stronger sorbent is nearest the pump. When used for personal sampling, the tube assembly shall be mounted vertically in the breathing zone. When used for fixed location sampling, a suitable sampling site is chosen.

Turn the pump on and adjust the flowrate so that the recommended sample air volume is taken in the available time. The recommended sample air volume for VOCs covered by this Standard is between 1 L and 10 L. If the total sample is likely to exceed 1 mg (1 mg on each tube), the sample air volume shall be reduced accordingly, or overload may occur.

Remarks 1 Sampling efficiency is 100 % (quantitative), provided that the sorptive capacity of the sorbents is not exceeded. If this capacity is exceeded, breakthrough of vapour from the tube assembly will occur. The breakthrough volume may be measured by sampling from a standard atmosphere, whilst monitoring the effluent air with a flame ionization or equivalent detector (a suitable method is described in Annex A). Alternatively, instead of determining the breakthrough volume directly, the mathematically related retention volume may be determined. The retention volume is determined chromatographically at elevated temperatures and subsequent extrapolation to room temperature. A suitable method is described in Annex B.

In the case where less hydrophobic strong sorbent is used for the measurement of the air which is 65 % or over of the relative humidity, a dehumidification tube filled with magnesium perchlorate should be mounted in front of the sorbent tube.

The breakthrough volume of porous polymers vary with ambient air temperature, reducing by a factor of about two for each 10 °C rise in temperature. It also varies with sampling flowrate, being reduced substantially at flowrates below 5 ml/min or above 500 ml/min. The breakthrough volumes of carbon molecular sieves are less affected by temperature and flowrate, but are substantially reduced at high concentrations of volatile organic vapour or high relative humidity. To allow a suitable margin of safety, a safe sampling volume (SSV) is defined such that it is a volume of not more than 70 % of the 5 % breakthrough volume (see **A.1.1** in Annex A) or 50 % of the retention volume (see in **B.1** in Annex B). Tables 1 to 6 give typical values for retention volumes and safe sampling volumes. These values have been determined by the chromatographic method (Annex B).

Remarks 2 The safe sampling volumes in tables 1 to 6 have been determined by the chromatographic method (Annex B). Measurements (**A.3**) by the direct method (Annex A) indicate that the chromatographic method is a reliable indication of the true breakthrough capacity except under conditions of high concentrations or very high humidity. These measurements (**A.3**) indicate that breakthrough volumes at high humidity (80 %) are about a factor of two lower for porous polymers and a factor of ten lower for carbonaceous sorbents than the low humidity value. If high concentrations [300 mg/m³ (100 ppm)] are also anticipated, the breakthrough volumes for carbonaceous sorbents should be further reduced by a factor of two.

If safe sampling volumes for compounds are estimated which are not listed in table 1, this estimation is only possible for such compounds which are situated between the two listed compounds of homologues of a chemical group. In all other cases the safe sampling volumes shall be tested experimentally with appropriate trials.

Note and record the times, temperature, or register reading if appropriate and the barometric pressure when the pump was turned on. At the end of the sampling period, note and record the flowrate or register reading, turn the pump off, and it is recommended to note and record the time, temperature and barometric pressure.

Disconnect the tube assembly and seal both ends of each tube with compression seals. Tighten these seals securely. The tubes should be uniquely labelled. Solvent-containing paints and markers or adhesive labels should not be used to label the tubes.

If samples are not to be analysed within 8 h, place them in a clean, uncoated, refrigerated sealed metal or glass container. If possible the sample should be refrigerated during transportation. Record air temperature and barometric pressure periodically during sampling if it is desired to express concentrations reduced to specific conditions (11.1).

Field blanks (travel blanks) should be prepared by using tubes identical to those used for sampling and subjecting them to the same handling procedure as the tubes except for the actual period of sampling. Label these as blanks.

Remarks 3 Since this method uses thermal desorption, unless the thermal desorption apparatus has the facility to retrap the sample after analysis, there will generally only be one opportunity to analyse the sample. If the sample is important and the chance of overload and/or sample breakthrough is a possibility, a second sample at a lower flowrate should be taken.

10 Procedure

10.1 Safety precautions This Standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this Standard to establish appropriate health and safety practices and determine the applicability of regulatory limitations prior to use.

10.2 Desorption and analysis Place the sorbent tube in a compatible thermal desorption apparatus. Purge the air from the tube to avoid chromatographic artefacts arising from the thermal oxidation of the sorbent or gas chromatographic stationary phase. Then heat the tube to displace the organic vapours which are passed to the gas chromatograph by means of a carrier gas stream. The gas flow direction at this stage should be the reverse of that used during sampling, i.e. the marked end of the tube should be nearest the gas chromatograph column inlet. Typically the gas flowrate through the tube should be in the order of 30 ml/min to 50 ml/min for optimum desorption efficiency.

For the initial air purge, it is usually necessary to use 10 times as much as the tube volume (i.e. 20 ml to 30 ml) of inert gas to completely displace the volume of air (2 ml to 3 ml) in the tube. However, if strongly hydrophilic sorbent are needed, it is necessary to employ a larger purge to reduce sorbed air and water to prevent ice formation blocking the cold trap. During the purge period, care should be taken to minimize heating of the tube.

The desorbed sample occupies a volume of several ml of gas, so that pre-concentration is essential prior to capillary gas chromatograph analysis. This can be achieved using a small, cooled, secondary sorbent trap, which can be desorbed sufficiently rapidly at low flowrates (<5 ml/min) to minimize band-broadening and produce capillary-compatible peaks. Alternatively, an empty secondary trap, or one containing an inert material such as glass beads, can be used to pre-concentrate the sample, but such traps typically require cooling to below -100°C . Alternatively, the desorbed sample can be passed directly to the gas chromatograph (single-stage desorption), where it shall be refocused. This typically requires a high phase-ratio column (e.g. 5 μm film thickness, 0.2 mm to 0.32 mm internal diameter) and a sub-ambient starting temperature.

If a secondary sorbent cold trap is not available and if sub-zero capillary cryofocusing temperatures are used to preconcentrate the analytes water shall be completely eliminated from the sample tube prior to desorption in order to prevent ice formation blocking the capillary tubing and stopping the thermal desorption process.

Remarks : If a secondary cold trap is not available and optimum sample tube desorption flowrates of 30 ml/min to 50 ml/min are used, a minimum split ratio of 30:1 to 50:1 will typically be required for operation with high-resolution capillary columns. Single-stage thermal desorption may thus limit sensitivity.

Desorption conditions should be chosen such that desorption from the sample tube is complete, and no sample loss occurs in the secondary trap, if used. Typical parameters are as follows:

Cold trap sorbent	typically same as tubes, 40 mg to 100 mg, if used
Cold trap temperature	$+20^{\circ}\text{C}$ to -180°C , depending on type of cold trap
Desorption temperature	250°C to 325°C
Desorption time	5 min to 15 min
Desorption flowrate	30 ml/min to 50 ml/min
Carrier gas	helium or nitrogen (selected according to the cold level)
Split ratios	Split ratios between the sample tube and secondary cold trap and between the secondary cold trap and analytical column (if applicable) shall be selected dependent on expected atmospheric concentration. (See guidance from respective manufacturers of the thermal desorption apparatus.)