

Status: Currently Official on 16-Feb-2025
 Official Date: Official as of 01-May-2021
 Document Type: USP Monographs
 DocId: GUID-0A83D713-83B4-4194-AF49-12B7F426D345_3_en-US
 DOI: https://doi.org/10.31003/USPNF_M67440_03_01
 DOI Ref: hr1jw

© 2025 USPC
 Do not distribute

Potassium Chloride in Lactated Ringer's and Dextrose Injection

DEFINITION

Potassium Chloride in Lactated Ringer's and Dextrose Injection is a sterile solution of Calcium Chloride, Potassium Chloride, Sodium Chloride, and Sodium Lactate in Water for Injection. It contains, in each 100 mL, NLT 285.0 and NMT 315.0 mg of sodium (Na) [as sodium chloride (NaCl) and anhydrous sodium lactate ($C_3H_5NaO_3$)], NLT 4.90 and NMT 6.00 mg of calcium (Ca) [equivalent to NLT 18.0 and NMT 22.0 mg of calcium chloride ($CaCl_2 \cdot 2H_2O$)], and NLT 231.0 and NMT 261.0 mg of lactate ($C_3H_5O_3$) [equivalent to NLT 290.0 and NMT 330.0 mg of anhydrous sodium lactate ($C_3H_5NaO_3$)]. It contains NLT 95.0% and NMT 105.0% of the labeled amount of potassium chloride (KCl), NLT 90.0% and NMT 105.0% of the labeled amount of dextrose ($C_6H_{12}O_6 \cdot H_2O$), and NLT 90.0% and NMT 110.0% of the labeled amount of chloride (Cl) [as sodium chloride (NaCl), potassium chloride (KCl), and calcium chloride ($CaCl_2 \cdot 2H_2O$)]. It contains no antimicrobial agents.

IDENTIFICATION

• A.

Sample solution: Nominally 50 mg/mL of dextrose from a suitable volume of Injection in [water](#)

Analysis: Add a few drops of *Sample solution* to 5 mL of hot [alkaline cupric tartrate TS](#).

Acceptance criteria: A copious red precipitate of cuprous oxide is formed.

Change to read:

• **B. IDENTIFICATION TESTS—GENERAL (191), [Chemical Identification Tests, Chloride](#) and [Calcium](#):** Meets the requirements of [▲]the [▲](USP 1-May-2021) test [▲](USP 1-May-2021) under *Chloride* and test A under *Calcium*

• **C. SODIUM:** The sample imparts an intense yellow color to a nonluminous flame.

• **D. POTASSIUM:** The sample imparts a violet color to a nonluminous flame. Because the presence of small quantities of sodium masks the color, screen out the yellow color produced by sodium by viewing through a blue filter that blocks the emission at 589 nm (sodium), but is transparent to emission at 404 nm (potassium). [NOTE—Traditionally, cobalt glass has been used, but other suitable filters are commercially available.]

• **E.** The retention time of the lactate peak of the *Sample solution* corresponds to that of the *Standard solution*, as obtained in the Assay for *Lactate*.

ASSAY

• CALCIUM

[NOTE—Concentrations of the *Standard solution* and the *Sample solution* may be modified to fit the linear or working range of the atomic absorption spectrophotometer.]

Lanthanum chloride solution: 88.45 g/L of [lanthanum chloride](#) prepared as follows. Transfer a suitable quantity of lanthanum chloride to an appropriate volumetric flask. Add 50% of the final flask volume of [water](#). Carefully add 25% of the final flask volume of [hydrochloric acid](#). Mix, and allow to cool. Dilute with [water](#) to volume.

Calcium stock solution: 1 mg/mL of calcium prepared as follows. Transfer 499.5 mg of primary standard calcium carbonate to a 200-mL volumetric flask, and add 10 mL of [water](#). Carefully add 5 mL of [diluted hydrochloric acid](#), and swirl to dissolve the calcium carbonate. Dilute with water to volume.

Standard solutions: 0.010, 0.015, and 0.020 mg/mL of calcium prepared as follows. To three separate 100-mL volumetric flasks, each containing 5.0 mL of *Lanthanum chloride solution*, add 1.0, 1.5, and 2.0 mL of *Calcium stock solution*, respectively. Dilute the contents of each flask with [water](#) to volume.

Sample solution: Transfer 30.0 mL of Injection to a 100-mL volumetric flask containing 5.0 mL of *Lanthanum chloride solution*. Dilute with [water](#) to volume.

Blank: Transfer 5.0 mL of *Lanthanum chloride solution* to a 100-mL volumetric flask and dilute with [water](#) to volume.

Instrumental conditions

(See [Atomic Absorption Spectroscopy \(852\)](#).)

Mode: Atomic absorption spectrophotometry

Analytical wavelength: Calcium emission line at 422.7 nm

Lamp: Calcium hollow-cathode

Flame: Air–acetylene

Analysis

Samples: *Standard solutions, Sample solution, and Blank*

Plot the absorbances of the *Standard solutions* versus the concentration, in mg/mL, of calcium, and draw the straight line best fitting the three plotted points. From the graph so obtained, determine the concentration (*C*), in mg/mL, of calcium in the *Sample solution*.

Calculate the quantity (mg) of calcium (Ca) in each 100 mL of Injection taken:

$$\text{Result} = C \times D \times F$$

C = concentration of calcium in the *Sample solution*, as determined from the graph (mg/mL)

D = dilution factor for the *Sample solution*, 3.3

F = conversion factor for each 100 mL of Injection, 100 mL

Acceptance criteria: 4.90–6.00 mg of calcium (Ca) in each 100 mL

• CHLORIDE

Sample solution: Transfer 10 mL of Injection into a conical flask. Add 10 mL of [glacial acetic acid](#), 75 mL of [methanol](#), and 3 drops of [eosin Y TS](#).

Titrimetric system

Mode: Direct titration

Titrant: [0.1 N silver nitrate VS](#)

Endpoint detection: Visual

Analysis

Sample: *Sample solution*

Titrate, with shaking, with *Titrant* to a pink endpoint.

Calculate the percentage of the labeled amount of chloride (Cl) in the portion of Injection taken:

$$\text{Result} = V \times N \times (F/W) \times 100$$

V = volume of *Titrant* consumed by the *Sample solution* (mL)

N = actual normality of the *Titrant* (mEq/mL)

F = equivalency factor, 35.45 mg/mEq

W = nominal amount of chloride in the *Sample solution* (mg)

Acceptance criteria: 90.0%–110.0%

• DEXTROSE

Sample solution: Transfer a volume of Injection containing 2–5 g of dextrose to a 100-mL volumetric flask. Add 0.2 mL of [6 N ammonium hydroxide](#), and dilute with [water](#) to volume.

Analysis

Sample: *Sample solution*

Determine the angular rotation in a suitable polarimeter tube (see [Optical Rotation \(781\)](#)).

Calculate the percentage of the labeled amount of dextrose ($\text{C}_6\text{H}_{12}\text{O}_6 \cdot \text{H}_2\text{O}$) in the portion of Injection taken:

$$\text{Result} = [(100 \times a)/(l \times \alpha)] \times (1/C_U) \times (M_{r1}/M_{r2}) \times 100$$

a = observed angular rotation of the *Sample solution* (°)

l = length of the polarimeter tube (dm)

α = midpoint of the specific rotation range for anhydrous dextrose, 52.9°

C_U = nominal concentration of dextrose in the *Sample solution* (g/100 mL)

M_{r1} = molecular weight of dextrose monohydrate, 198.17

M_{r2} = molecular weight of anhydrous dextrose, 180.16

Acceptance criteria: 90.0%–105.0%