



Association of Textile, Apparel & Materials Professionals

AATCC takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this test method. Users of this test method are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.

This test method is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reaffirmed or withdrawn. Your comments are invited either for revision of this test method or for additional procedures and should be addressed to the AATCC Technical Center. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing, you should make your views known to Diana A. Wyman, Technical Director, at the address shown below.

This test method is copyrighted by AATCC, PO Box 12215, Research Triangle Park, NC 27709 USA. Individual reprints (single or multiple copies) of this method may be obtained by contacting AATCC at the above address or tel: 919-549-3526; fax: 919-549-8933, or e-mail: orders@aatcc.org.

AATCC License Agreement: This test method is copyrighted by the American Association of Textile Chemists and Colorists (AATCC), PO Box 12215, 1 Davis Drive, Research Triangle Park, NC USA. All rights reserved.

AATCC grants you a license as follows: The right to download one electronic file of this AATCC test method for temporary storage on one computer for purposes of viewing, and/or printing one copy of the AATCC test method for individual use. Neither the electronic file nor the hard copy print may be reproduced in any way. In addition, the electronic file may not be distributed elsewhere over computer networks or otherwise. The hard copy print may only be distributed to other employees for their internal use within your organization.

This test method is not for resale.

Chelation Value of Polyaminopolycarboxylic Acids and Their Salts: Copper PAN Method

Developed in 1987 by AATCC Committee RR90; reaffirmed 1988 (with title change), 2002, 2012; editorially revised and reaffirmed 1989, 1997, 2007; revised 1992, 2018 (with title change); editorially revised 2010.

Foreword

This test method is an alternate to AATCC Test Method (TM) 149, Chelation Value of Aminopolycarboxylic Acids and Their Salts: Calcium Oxalate Method. Titration with a calcium salt using either oxalate or a dye indicator, as in TM149, yields a value that includes partially-substituted EDTA, HEDTA, DTPA, nitrilotriacetic acid (NTA), iminodiacetic acid (IDA), glycolate and other weak chelating agents found in some commercial products. This alternative procedure excludes the effect of the partially-substituted products, providing a more stringent value.

1. Purpose and Scope

1.1 The active ingredient content of ethylenediaminetetraacetic acid (EDTA), N-hydroxyethylethylenediaminetriacetic acid (HEDTA) and diethylenetriaminepentaacetic acid (DTPA) and their salts is usually expressed by the calcium chelation value (CaCV). This value indicates the amount of calcium (as calcium carbonate) which will be chelated by a known weight of chelating agent. This method excludes the effect of partially-substituted products.

2. Principle

2.1 Calcium chelation value is determined experimentally by titrating a weighed specimen of chelating agent with a copper nitrate solution of known concentration in the presence of PAN (1-(2-pyridylazo)-2-naphthol). Initially, the color of the solution is yellow-green due to free indicator. When all of the chelating agent has reacted with copper nitrate, the color of the solution changes to a permanent purple which is the end point.

2.2 The CaCV indicates the amount of calcium (as calcium carbonate) that will be sequestered by a known weight of chelating agent and is expressed as milligrams of calcium carbonate (CaCO_3) per gram of chelating agent (mg calcium carbonate/g chelating agent).

3. Terminology

3.1 **chelating agent**, n.—in *textile chemistry*, a chemical capable of deactivating metal ions by formation of a water-soluble complex. SYN: *sequestering agent*.

4. Safety Precautions

NOTE: These safety precautions are for information purposes only. The precautions are ancillary to the testing procedures and are not intended to be all inclusive. It is the user's responsibility to use safe and proper techniques in handling materials in this test method. Manufacturers MUST be consulted for specific details such as material safety data sheets and other manufacturer's recommendations. All OSHA standards and rules must also be consulted and followed.

4.1 Good laboratory practices should be followed. Wear safety glasses in all laboratory areas.

4.2 All chemicals should be handled with care.

4.3 In preparing, dispensing and handling glacial acetic acid and sodium hydroxide use chemical goggles or face shield, impervious gloves and an impervious apron. Concentrated acids and bases should be handled only in an adequately ventilated laboratory hood. CAUTION: Always add acid to water.

4.4 Cupric nitrate 2.5 hydrate is corrosive to the eye and skin and is harmful if inhaled. It is an oxidizing material and may react with organic materials. It should be handled only in an adequately ventilated hood.

4.5 Methanol is harmful if inhaled or swallowed and is a flammable liquid. It should be stored in the laboratory only in small containers away from heat, open flames and sparks. This chemical should not be used near an open flame and should be handled only in an adequately ventilated hood.

4.6 An eyewash/safety shower should be located nearby and a self-contained breathing apparatus should be readily available for emergency use.

4.7 Exposure to chemicals used in this procedure must be controlled at or below levels set by governmental authorities (e.g., Occupational Safety and Health Administration's [OSHA] permissible exposure limits [PEL] as found in 29 CFR 1910.1000; see www.osha.gov for latest

version). In addition, the American Conference of Governmental Industrial Hygienists (ACGIH) Threshold Limit Values (TLVs) comprised of time weighted averages (TLV-TWA), short term exposure limits (TLV-STEL) and ceiling limits (TLV-C) are recommended as a general guide for air contaminant exposure which should be met (see 14.1).

5. Reagents

- 5.1 Acetic acid, 1.0N (CH_3COOH)
- 5.2 Cupric nitrate 2.5 hydrate ($\text{Cu}(\text{NO}_3)_2 \cdot 2.5\text{H}_2\text{O}$)
- 5.3 Ethylene diaminetetraacetic acid (free acid) (EDTA) ($\text{C}_{10}\text{H}_{16}\text{O}_8\text{N}_2$)
- 5.4 Methanol
- 5.5 1-(2-Pyridylazo)-2-naphthol, (PAN) ($\text{C}_{15}\text{H}_{11}\text{ON}_3$)
- 5.6 Sodium acetate trihydrate ($\text{CH}_3\text{COONa} \cdot 3\text{H}_2\text{O}$)
- 5.7 Sodium hydroxide, 1.0N (NaOH)

6. Sampling

6.1 Conduct the test in triplicate; i.e., analyze 3 test specimens from a sample.

7. Specimens

7.1 Do not use aluminum or metallic weighing pans.

7.2 For analysis of solid EDTA, HEDTA, and DTPA (either free acid or salt forms) weigh a 0.24-0.26 g specimen of the dried chelating agent to an accuracy of 0.01 g.

7.3 For analysis of commercial solution, of salts of the chelating agents, weigh a 0.49-0.51 g specimen of the solution to an accuracy of 0.01 g.

8. Conditioning

8.1 If solid chelating agent is being evaluated, dry a 2-g sample at the appropriate temperature for at least 2 h and cool in a desiccator before weighing specimens.

8.2 Dry free acid forms of chelating agents at 120°C.

8.3 Dry salt forms of chelating agents at 80°C.

9. Preparation of Reagents

9.1 Cupric Nitrate Solution: Dissolve 23.30 g of cupric nitrate 2.5 hydrate and dilute to 1.000 L in a volumetric flask.