



Designation: E 961 – 97

AMERICAN SOCIETY FOR TESTING AND MATERIALS  
100 Barr Harbor Dr., West Conshohocken, PA 19428  
Reprinted from the Annual Book of ASTM Standards. Copyright ASTM

## Standard Specification for Blood Sedimentation Tube, Wintrobe, Glass, Reusable<sup>1</sup>

This standard is issued under the fixed designation E 961; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ε) indicates an editorial change since the last revision or reapproval.

### 1. Scope

1.1 This specification covers reusable blood sedimentation tubes suitable for determining sedimentation rates and the volume of packed red blood cells.

1.2 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.*

### 2. Referenced Documents

#### 2.1 ASTM Standards:

- E 438 Specification for Glasses in Laboratory Apparatus<sup>2</sup>
- E 671 Specification for Maximum Permissible Thermal Residual Stress in Annealed Glass Laboratory Apparatus<sup>2</sup>
- E 920 Specification for Commercially Packaged Laboratory Apparatus<sup>2</sup>
- E 921 Specification for Export Packaged Laboratory Apparatus<sup>2</sup>
- E 1133 Practice for Performance Testing of Packaged Laboratory Apparatus for United States Government Procurements<sup>2</sup>
- E 1157 Specification for Sampling and Testing of Reusable Laboratory Glassware<sup>2</sup>

### 3. Terminology

#### 3.1 Definitions of Terms Specific to This Standard:

- 3.1.1 *reusable*—Capable of being used again.
- 3.1.2 *Wintrobe*—The surname of the individual responsible for the design of the Wintrobe tube and the method of use.<sup>3</sup>

### 4. Classification

4.1 This specification covers a tube that is intended to be used until it is no longer considered a functional device for the purpose intended.

### 5. Materials

5.1 *Glass*—The tubes made to this specification shall be

fabricated from borosilicate glass, Type I, Class B, or soda-lime glass, Type II, in accordance with Specification E 438.

### 6. Dimensions and Graduations

6.1 *Dimensions*—The tube shall be made of tubing with an outside diameter (O.D.) of 7.0 to 8.0 mm with an inside diameter (I.D.) of 2.9 to 3.3 mm. The uniformity of the bore shall be  $\pm 0.1$  mm throughout the tube. The tube shall be 110 to 117 mm long and have a graduated scale of  $105 \pm 0.25$  mm from the inside bottom of the tube. The tube shall be legibly marked with the manufacturer's or vendor's name or mark and possess a frosted area for marking purposes.

6.2 *Graduation Scale*—The tube shall be graduated  $105 \pm 0.25$  mm in 1-mm divisions and numbered every 1 cm with two sets of numerals. One set of graduation numerals shall be from 0 to 9 cm down the left side of the graduation scale and the other set of graduation numerals shall be from 1 to 10 cm up the right side of the graduation scale as shown in Fig. 1.

6.3 *Graduation Lines*—The graduation lines shall be of uniform width with a maximum width of 0.25 mm for acid-etched and filled-in lines and 0.4 mm for amber-stain, fired-in lines. They shall be at right angles to the vertical axis of the tube with a maximum of 0.2 mm tolerance between two adjacent lines. Every tenth (numbered) line shall be a minimum of 6 mm long. The medium (every fifth) line shall be a minimum of 4 mm long. The short (intermediate) lines shall be a minimum of 2.5 mm long.

### 7. Workmanship, Finish, and Appearance

7.1 *Design*—The tube shall be made of thick-walled glass tubing. The tube shall have a cylindrical body terminating in a round bottom on the outside and a flat bottom on the inside. The tube shall have a uniform bore and a beaded flared top.

7.2 *Workmanship*—The tube shall be free of physical defects which will impair its use. The tube shall be well annealed with no ring strain when tested as specified in 8.1. Longitudinal strain, if present, shall be faint and highly diffused. The tube shall not chip, crack, or break during centrifugation when tested as specified in 8.3.

7.3 *Marking Permanency*—Inscriptions, graduation lines and numerals shall be either etched and filled with a permanent pigment, or an amber stain fired into the glass tube. The color depth of the markings on the tube shall be adequate to permit routine functional use of the tube without creating a difficulty in setting a meniscus or reading the separation of blood cells


<sup>1</sup> This specification is under the jurisdiction of ASTM Committee E-41 on Laboratory Apparatus and is the direct responsibility of Subcommittee E 41.01 on Glass Apparatus.

Current edition approved Nov. 10, 1997. Published February 1998. Originally published as E 961 – 83. Last previous edition E 961 – 93.

<sup>2</sup> *Annual Book of ASTM Standards*, Vol 14.02.

<sup>3</sup> Wintrobe, Maxwell M., "Laboratory Evaluation of Erythrocytes," *Clinical Hematology*, Seventh Ed., 1974, pp. 109–134.

**ASTM E 961**



| FEATURE          |   | DIMENSIONS IN mm         |
|------------------|---|--------------------------|
| LENGTH           | A | 110 – 117                |
| OUTSIDE DIAMETER | B | 7.0 – 8.0                |
| INSIDE DIAMETER  | C | 2.9 – 3.3                |
| SCALE LENGTH     | D | 105 ± .25                |
| FROSTED          | E | APPROX. 6.0 O.D. ON BACK |

**FIG. 1 Wintrobe Tube**

from plasma. The permanency of the markings shall meet the requirements of the test described in 8.2.

## 8. Test Methods

**8.1 Strain Examination**—Test tubes in accordance with the polarimeter measurement procedure as outlined in Specification E 671.

**8.2 Pigmentation and Amber Stain Permanency Test**—Freshly prepare a chromic acid cleaning solution by combining 200 g of sodium dichromate ( $\text{Na}_2\text{Cr}_2\text{O}_7 \times 2\text{H}_2\text{O}$ ), 1 000 mL of distilled water and 1 500 mL of sulfuric acid ( $\text{H}_2\text{SO}_4$ , ACS Reagent 95 to 98 %). Immerse the tube in the chromic acid solution. Let stand at room temperature (20 to 25°C) for 15 min. Remove tube from the solution and thoroughly rinse in tap water followed by distilled water. Dry the tube by rubbing vigorously, 5 to 10 strokes, with a laboratory cloth or tissue. This appearance of the markings should be the same as before the test, when judged by eye under normal room lighting.

**8.3 Resistance to Centrifugal Force Test**—Fill the tube with water to the top graduation line and place in a centrifuge. The

speed and dimensions of centrifuge head shall be such that the inside bottom of the tube is subjected to a relative centrifugal force of (RCF) not less than 2 500 gravities. Calculate the relative centrifugal force as follows:

$$\begin{aligned} RCF &= \text{relative centrifugal force in gravities,} \\ RCF &= 00001118 \times r \times N^2 \text{ gravities,} \\ r &= \text{rotating radius to inside or outside of tube in centimetre,} \\ &\text{and} \\ N &= \text{rotating speed in revolutions per minute.} \end{aligned}$$

8.4 For additional sampling and testing data, see Specification E 1157.

## 9. Packaging

9.1 For packaging, select from either Specification E 920, Specification E 921, or Practice E 1133.

## 10. Keywords

10.1 blood; glass; reusable; sedimentation rate; wintrobe

*The American Society for Testing and Materials takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this standard. Users of this standard are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.*

*This standard is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reapproved or withdrawn. Your comments are invited either for revision of this standard or for additional standards and should be addressed to ASTM Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing you should make your views known to the ASTM Committee on Standards, 100 Barr Harbor Drive, West Conshohocken, PA 19428.*