

Status: Currently Official on 17-Feb-2025
Official Date: Official as of 01-May-2018
Document Type: USP Monographs
DocId: GUID-86A09EF6-8317-42F2-853F-7BEC550862DF_3_en-US
DOI: https://doi.org/10.31003/USPNF_M6730_03_01
DOI Ref: 3hq3g

© 2025 USPC
Do not distribute

Azathioprine Sodium for Injection

DEFINITION

Azathioprine Sodium for Injection is a sterile solid prepared by the freeze-drying of an aqueous solution of Azathioprine and Sodium Hydroxide. It contains NLT 93.0% and NMT 107.0% of the labeled amount of azathioprine ($C_9H_7N_7O_2S$).

IDENTIFICATION

- A.** The principal spot from the *Sample solution* shows the same R_F value as that obtained from *Standard solution A* in the test for *Limit of Mercaptopurine*.
- B.** The retention time of the major peak of the *Sample solution* corresponds to that of the *Standard solution*, as obtained in the Assay.

ASSAY

PROCEDURE

Solution A: 1.6 g/L of [sodium 1-heptanesulfonate](#) in [water](#)

Mobile phase: [Methanol](#) and *Solution A* (30:70). Adjust with [1 N hydrochloric acid](#) to a pH of 3.5 ± 0.1 .

Standard stock solution: 0.5 mg/mL of [USP Azathioprine RS](#) prepared as follows. Transfer [USP Azathioprine RS](#) to a suitable volumetric flask. Add 25% of the flask volume of [methanol](#) and 1% of [ammonium hydroxide](#) to the flask, swirl, and sonicate for 2 min or until dissolved. Dilute with [methanol](#) to volume.

Standard solution: 0.1 mg/mL of [USP Azathioprine RS](#) in [water](#) from the *Standard stock solution*

Sample solution: Nominally equivalent to 0.1 mg/mL of azathioprine in [water](#) from Azathioprine Sodium for Injection

Chromatographic system

(See [Chromatography \(621\), System Suitability](#).)

Mode: LC

Detector: UV 254 nm

Column: 3.9-mm $\times$ 30-cm; 10- μ m packing [L1](#)

Flow rate: 1.8 mL/min

Injection volume: 10 μ L

System suitability

Sample: *Standard solution*

Suitability requirements

Tailing factor: NMT 2

Relative standard deviation: NMT 1.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the percentage of the labeled amount of azathioprine ($C_9H_7N_7O_2S$) in the portion of Azathioprine Sodium for Injection taken:

$$\text{Result} = (r_U/r_S) \times (C_S/C_U) \times 100$$

r_U = peak response of azathioprine from the *Sample solution*

r_S = peak response of azathioprine from the *Standard solution*

C_S = concentration of [USP Azathioprine RS](#) in the *Standard solution* (mg/mL)

C_U = nominal concentration of azathioprine in the *Sample solution* (mg/mL)

Acceptance criteria: 93.0%–107.0%

PERFORMANCE TESTS

- UNIFORMITY OF DOSAGE UNITS (905):** Meets the requirements

IMPURITIES

LIMIT OF MERCATOPURINE

Standard solution A: 10 mg/mL of [USP Azathioprine RS](#) in [dimethylformamide](#)

Standard solution B: 100 μ g/mL of [USP Mercaptopurine RS](#) in [dimethylformamide](#)

Sample solution: 10 mg/mL of Azathioprine Sodium for Injection in [dimethylformamide](#)

Chromatographic system

(See [Chromatography \(621\), Thin-Layer Chromatography](#).)

Mode: TLC

Adsorbent: 0.25-mm layer of microcrystalline cellulose

Application volume: 5 µL for *Standard solution A* and the *Sample solution*, and 15 µL for *Standard solution B*

Developing solvent system: [Butyl alcohol](#) saturated with 5 N [ammonium hydroxide](#)

Analysis

Samples: *Standard solution A*, *Standard solution B*, and *Sample solution*

Proceed as directed in the chapter. Apply the *Samples* at points 2 cm from the bottom edge of a TLC plate. Allow the spots to dry, and develop the chromatogram in a suitable chamber until the solvent front has moved 15 cm from the point of application. Remove the plate, air-dry, and locate the spots by viewing under short- and long-wavelength UV light.

Acceptance criteria: 3.0%; any spot from the *Sample solution*, other than the principal spot, is not more intense than the spot from *Standard solution B*.

SPECIFIC TESTS

• **pH (791).**

Sample solution: The contents of one container dissolved in 10 mL of [water](#)

Acceptance criteria: 9.8–11.0

• **BACTERIAL ENDOTOXINS TEST (85):** It contains NMT 1.0 USP Endotoxin Unit/mg of azathioprine.

• **WATER DETERMINATION (921), Method I:** NMT 7.0%, when the *Test Preparation* is prepared as directed for a hygroscopic specimen

• **OTHER REQUIREMENTS:** It meets the requirements in [Injections and Implanted Drug Products \(1\)](#).

ADDITIONAL REQUIREMENTS

• **PACKAGING AND STORAGE:** Preserve as described in [Packaging and Storage Requirements \(659\)](#), [Injection Packaging](#), [Packaging for constitution](#), at controlled room temperature.

• **USP REFERENCE STANDARDS (11).**

[USP Azathioprine RS](#)

[USP Mercaptopurine RS](#)

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
AZATHIOPRINE SODIUM FOR INJECTION	Documentary Standards Support	SM32020 Small Molecules 3
REFERENCE STANDARD SUPPORT	RS Technical Services RSTECH@usp.org	SM32020 Small Molecules 3

Chromatographic Database Information: [Chromatographic Database](#)

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. PF 41(6)

Current DocID: GUID-86A09EF6-8317-42F2-853F-7BEC550862DF_3_en-US

Previous DocID: GUID-86A09EF6-8317-42F2-853F-7BEC550862DF_1_en-US

DOI: https://doi.org/10.31003/USPNE_M6730_03_01

DOI ref: [3hq3g](#)