



BSI Standards Publication

Soaps — Determination of chloride content — Potentiometric method

National foreword

This British Standard is the UK implementation of [ISO 4323:2018](#). It supersedes [BS 1715-2.7:1989](#), which is withdrawn.

The UK participation in its preparation was entrusted to Technical Committee CII/34, Methods of test for surface active agents.

A list of organizations represented on this committee can be obtained on request to its secretary.

This publication does not purport to include all the necessary provisions of a contract. Users are responsible for its correct application.

© The British Standards Institution 2018
Published by BSI Standards Limited 2018

ISBN 978 0 580 92742 3

ICS 71.100.40

Compliance with a British Standard cannot confer immunity from legal obligations.

This British Standard was published under the authority of the Standards Policy and Strategy Committee on 31 December 2018.

Amendments/corrigenda issued since publication

Date	Text affected
------	---------------

INTERNATIONAL STANDARD

ISO
4323

Second edition
2018-12-13

Soaps — Determination of chloride content — Potentiometric method

Savons — Dosage des chlorures — Méthode potentiométrique



Reference number
ISO 4323:2018(E)
Documents
Document Sharing Hub
© ISO 2018



COPYRIGHT PROTECTED DOCUMENT

© ISO 2018, Published in Switzerland

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized otherwise in any form or by any means, electronic or mechanical, including photocopying, or posting on the internet or an intranet, without prior written permission. Permission can be requested from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
Ch. de Blandonnet 8 • CP 401
CH-1214 Vernier, Geneva, Switzerland
Tel. +41 22 749 01 11
Fax +41 22 749 09 47
copyright@iso.org
www.iso.org

Contents

Page

Foreword	iv
1 Scope	1
2 Normative references	1
3 Terms and definitions	1
4 Principle	1
5 Reagents	1
6 Apparatus	2
7 Sampling	2
8 Procedure	2
8.1 Measurement temperature.....	2
8.2 Calibration of the silver nitrate solution.....	3
8.2.1 Titration.....	3
8.2.2 Calculation of normality of solution.....	4
8.3 Determination.....	4
8.3.1 Test portion.....	4
8.3.2 Titration.....	4
8.3.3 Blank test.....	4
9 Expression of results	5
10 Test report	5
Annex A (informative) Example of how to calculate V_{EQ}	6