

Status: Currently Official on 15-Feb-2025
Official Date: Official as of 01-May-2023
Document Type: USP Monographs
DocId: GUID-5191DADB-384D-4AE7-BB40-DD48273667D3_6_en-US
DOI: https://doi.org/10.31003/USPNF_M47150_06_01
DOI Ref: an99b

© 2025 USPC
Do not distribute

Change to read:

Manganese Sulfate Injection

▲ (USP 1-May-2023)

DEFINITION

Manganese Sulfate Injection is a sterile solution of Manganese Sulfate in Water for Injection. It contains NLT 90.0% and NMT 110.0% of the labeled amount of manganese (Mn).

IDENTIFICATION

• **A.** The *Sample solution*, prepared as directed in the Assay, exhibits an absorption maximum at about 279 nm when tested as directed for in the Assay.

ASSAY

• PROCEDURE

Solution A: 0.9 mg/mL of [sodium chloride](#)

Standard stock solution A: 1 mg/mL of [manganese](#) prepared as follows. Transfer 1.000 g of [manganese](#) to a 1000-mL volumetric flask, dissolve in 20 mL of [nitric acid](#), and dilute with [0.1 N hydrochloric acid](#) to volume. Store in a polyethylene bottle.

Standard stock solution B: 20 µg/mL of manganese from *Standard stock solution A* in [water](#)

Standard solution A: 1.6 µg/mL of [manganese](#) prepared as follows. Transfer 4.0 mL of *Standard stock solution B* to a 50-mL volumetric flask containing 10 mL of *Solution A*, and dilute with [water](#) to volume.

Standard solution B: 2.0 µg/mL of [manganese](#) prepared as follows. Transfer 5.0 mL of *Standard stock solution B* to a 50-mL volumetric flask containing 10 mL of *Solution A*, and dilute with [water](#) to volume.

Standard solution C: 2.4 µg/mL of [manganese](#) prepared as follows. Transfer 6.0 mL of *Standard stock solution B* to a 50-mL volumetric flask containing 10 mL of *Solution A*, and dilute with [water](#) to volume.

Sample solution: Nominally 0.002 mg/mL of manganese from a volume of Injection in [water](#)

Instrumental conditions

(See [Atomic Absorption Spectroscopy \(852\)](#).)

Mode: Atomic absorption spectrophotometer

Analytical wavelength: Manganese emission line at 279 nm

Source: Manganese hollow-cathode lamp

Flame: Air–acetylene

Blank: 0.18 mg/mL of [sodium chloride](#)

Analysis

Samples: *Standard solution A*, *Standard solution B*, *Standard solution C*, *Sample solution*, and *Blank*

Plot the absorbances of the *Standard solutions* versus concentration, in µg/mL, of manganese, and draw the straight line best fitting the three plotted points. From the graph so obtained, determine the concentration, C_s , in µg/mL, of manganese in the *Sample solution*.

Calculate the percentage of the labeled amount of manganese in the portion of Injection taken:

$$\text{Result} = (C_s/C_U) \times 100$$

C_s = concentration of manganese in the *Sample solution*, determined from the calibration graph (µg/mL)

C_U = nominal concentration of manganese in the *Sample solution* (µg/mL)

Acceptance criteria: 90.0%–110.0%

SPECIFIC TESTS

Change to read:

- **BACTERIAL ENDOTOXINS TEST (85):** ▲Meets the requirements▲ (USP 1-May-2023)
- **pH (791):** 2.0–3.5
- **PARTICULATE MATTER IN INJECTIONS (788):** It meets the requirements for small-volume injections.
- **OTHER REQUIREMENTS:** It meets the requirements under [Injections and Implanted Drug Products \(1\)](#).

ADDITIONAL REQUIREMENTS

- **PACKAGING AND STORAGE:** Preserve in single-dose or in multiple-dose containers, preferably of Type I or Type II glass.
- **LABELING:** Label the Injection to indicate that it is to be diluted to the appropriate strength with Sterile Water for Injection or other suitable fluid prior to administration.

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
MANGANESE SULFATE INJECTION	Documentary Standards Support	SM32020 Small Molecules 3

Chromatographic Database Information: [Chromatographic Database](#)

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. 47(5)

Current DocID: GUID-5191DADB-384D-4AE7-BB40-DD48273667D3_6_en-US

DOI: https://doi.org/10.31003/USPNF_M47150_06_01

DOI ref: [an99b](#)

OFFICIAL