


Status: Currently Official on 16-Feb-2025
 Official Date: Official as of 01-Aug-2023
 Document Type: NF Monographs
 DocId: GUID-A86A8C0D-EAE1-4672-82F9-A9ED2FA56C0D_3_en-US
 DOI: https://doi.org/10.31003/USPNF_M15735_03_01
 DOI Ref: u1t0x


© 2025 USPC
 Do not distribute

Add the following:

▲Polyethylene Glycol 12 Cetostearyl Ether



or



$n = 12$ (average)

Alcohols, C16–C18, ethoxylated

CAS RN[®]: 68439-49-6.

DEFINITION

Polyethylene Glycol 12 Cetostearyl Ether is a mixture of ethers of mixed polyethylene glycols with linear fatty alcohols, mainly cetostearyl alcohol. It may contain some free polyethylene glycols and various amounts of free cetostearyl alcohol. The number of moles of ethylene oxide reacted per mole of cetostearyl alcohol is 12 (nominal value).

IDENTIFICATION

- **A. SPECTROSCOPIC IDENTIFICATION TESTS (197), Infrared Spectroscopy:** 197A
- **B. ▲HYDROXYL VALUE: ▲** (ERR 1-Aug-2023) It meets the requirements in *Specific Tests for [Fats and Fixed Oils \(401\)](#), Procedures, Hydroxyl Value*.
- **C. TEST FOR STEARYL ALCOHOL**

Standard solution: Dissolve 25 mg of [USP Stearyl Alcohol RS](#) in [methanol](#) and dilute with the same solvent to 25 mL.

Sample solution: Dissolve about 10.0 g of the substance to be examined in a mixture of 1 volume of water and 9 volumes of [methanol](#) and dilute with the same mixture of solvents to 75 mL. Add 60 mL of [hexane](#) and shake for 3 min. The formation of foam can be reduced by the addition of 2–3 drops of [alcohol](#). Filter the upper layer through [anhydrous sodium sulfate](#), wash the filter with 3 quantities, each with 10 mL of [hexane](#), and evaporate the combined filtrates to dryness. Dissolve 0.05 g of the residue in 10 mL of [methanol](#) (the solution may be opalescent).

Chromatographic system

(See [Chromatography \(621\)](#), *General Procedures, Thin-Layer Chromatography* and [Thin-Layer Chromatographic Identification Test \(201\)](#).)

Mode: TLC

Adsorbent: Chromatographic silica gel with an average particle size of 5 µm (HPTLC plate)¹

Developing solvent system: [Ethyl acetate](#)

Application volume: 20 µL

Temperature: Ambient

Developing distance: 15 cm

Spray reagent: Dissolve 0.5 g of *vanillin* in 50 mL of [alcohol](#), and dilute with [sulfuric acid](#) to 100 mL.

Analysis

Samples: Apply the *Standard solution* and *Sample solution*, and allow the spots to dry. Develop the chromatogram until the solvent front has moved about 15 cm. Remove the plate from the developing chamber, mark the solvent front, and allow the solvent to evaporate. Spray with *Spray reagent*; allow to dry in air; heat at about 130° for 15 min, and allow to cool in air.

Acceptance criteria: The chromatogram obtained with the *Sample solution* shows several spots. The retardation factor (R_f) of one of these spots corresponds to the R_f of the principal spot in the chromatogram obtained with the *Standard solution*.



- **D.**
Diluted hydrochloric acid: Dilute 20 g of [hydrochloric acid](#) with water to 100 mL.
Barium chloride solution: Prepare a 61-g/L solution of [barium chloride](#) in water.
Phosphomolybdic acid solution: Prepare a 100-g/L solution of [phosphomolybdic acid](#) in water.
Analysis: Dissolve or disperse the Sample in ▲5 mL▲ (ERR 1-Aug-2023) [alcohol](#). Add 2 mL of water, 10 mL of *Diluted hydrochloric acid*, 10 mL of *Barium chloride solution*, and 10 mL of *Phosphomolybdic acid solution*.
Sample: 0.1 g
Acceptance criteria: A precipitate is formed.

IMPURITIES

- **RESIDUE ON IGNITION (281):** NMT 0.2%
- **ETHYLENE OXIDE AND DIOXANE (228), Method I**
Acceptance criteria
For ethylene oxide: NMT 1 ppm
For dioxane: NMT 10 ppm

SPECIFIC TESTS

- **FATS AND FIXED OILS (401), Procedures, Acid Value, Method I**
Analysis: Accurately weigh 5.0 g of the sample. [NOTE—[Petroleum ether](#) with a 100°–120° boiling range can be used to replace [ether](#) in the test.]
Acceptance criteria: NMT 1.0
- **FATS AND FIXED OILS (401), Procedures, Hydroxyl Value:** 67–77
- **FATS AND FIXED OILS (401), Procedures, Saponification Value**
Analysis: Accurately weigh 10.0 g of the sample.
Acceptance criteria: NMT 3.0
- **FATS AND FIXED OILS (401), Procedures, Iodine Value, Method I**
Analysis: Accurately weigh 1.0 g of the sample.
Acceptance criteria: NMT 2.0
- **WATER DETERMINATION (921), Method I**
Analysis: Accurately weigh 2.0 g of the sample.
Acceptance criteria: NMT 3.0%
- **ALKALINITY**
Sample: 2.0 g
Bromothymol blue solution: Dissolve 50 mg of [bromothymol blue](#) indicator in a mixture of 4 mL of [0.02 N sodium hydroxide](#) and 20 mL of [alcohol](#), and dilute with water to 100 mL.
Analysis: Dissolve the sample in a hot mixture of 10 mL of [alcohol](#) and 10 mL of water. Add 0.1 mL of *Bromothymol blue solution*.
Acceptance criteria: NMT 0.5 mL of [0.1 N hydrochloric acid](#) is required to change the color of the indicator to yellow.

ADDITIONAL REQUIREMENTS

- **PACKAGING AND STORAGE:** Preserve in tight containers, at ambient temperature.
- **USP REFERENCE STANDARDS (11)**
[USP Polyethylene Glycol 12 Cetostearyl Ether RS](#)
[USP Stearyl Alcohol RS](#)▲ (NF 1-Aug-2023)

¹ A suitable grade is available as catalog #1.05641.0001 from [www.EMDMillipore.com](#).

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
POLYETHYLENE GLYCOL 12 CETOSTEARYL ETHER	Documentary Standards Support	CE2020 Complex Excipients  Documents Document Sharing Hub