

Status: Currently Official on 16-Feb-2025
 Official Date: Official as of 01-Aug-2023
 Document Type: USP Monographs
 DocId: GUID-5340F03C-9CCF-4D61-ADC7-A7A9C2AD7AC9_3_en-US
 DOI: https://doi.org/10.31003/USPNF_M72770_03_01
 DOI Ref: c9hdz

© 2025 USPC
 Do not distribute

Quinidine Sulfate Tablets

DEFINITION

Quinidine Sulfate Tablets contain amounts of quinidine sulfate and dihydroquinidine sulfate totaling NLT 90.0% and NMT 110.0% of the labeled amount of quinidine sulfate, calculated as quinidine sulfate $[(C_{20}H_{24}N_2O_2)_2 \cdot H_2SO_4 \cdot 2H_2O]$.

IDENTIFICATION

• A.

Analysis: Shake the equivalent to 100 mg of quinidine sulfate from powdered Tablets, with 10 mL of dilute sulfuric acid (1 in 350), and filter.

Acceptance criteria: An appropriate dilution of the filtrate exhibits a vivid blue fluorescence. On the addition of a few drops of hydrochloric acid, the fluorescence disappears.

• B. The R_f value of the principal spot of the *Sample solution* corresponds to that of the *Standard solution*, as obtained in *Organic Impurities*.

• C. [IDENTIFICATION TESTS—GENERAL, Sulfate\(191\)](#).

Analysis: Shake the equivalent to 100 mg of quinidine sulfate from powdered Tablets, with 10 mL of dilute hydrochloric acid (1 in 100), and filter.

Acceptance criteria: The filtrate meets the requirements.

ASSAY

• PROCEDURE

Solution A: Add 35.0 mL of methanesulfonic acid to 20.0 mL of glacial acetic acid, and dilute with water to 500 mL.

Solution B: Dilute 10.0 mL of diethylamine with water to 100 mL.

Mobile phase: Acetonitrile, *Solution B*, *Solution A*, and water (100:20:20:860). Adjust with diethylamine to a pH of 2.6.

System suitability solution: 0.2 mg/mL each of [USP Quinidine Sulfate RS](#) and dihydroquinidine hydrochloride prepared as follows. Dissolve a suitable quantity of each [USP Quinidine Sulfate RS](#) and dihydroquinidine hydrochloride in about 10% of total volume of methanol, and dilute with *Mobile phase* to volume.

Standard solution: 0.2 mg/mL of [USP Quinidine Sulfate RS](#) in *Mobile phase*

Sample stock solution: Weigh and finely powder NLT 20 Tablets. Transfer a portion of the powder, nominally equivalent to 160 mg of quinidine sulfate, to a 100-mL volumetric flask. Add 80 mL of methanol, and shake the flask by mechanical means for 30 min. Dilute with methanol to volume, and filter, discarding the first 10 mL of filtrate.

Sample solution: Nominally equivalent to 0.192 mg/mL of quinidine sulfate in *Mobile phase* from the *Sample stock solution*

Chromatographic system

(See [Chromatography \(621\)](#), [System Suitability](#).)

Mode: LC

Detector: UV 235 nm

Column: 3.9-mm × 30-cm; packing L1

Injection volume: 50 µL

System suitability

Sample: *System suitability solution*

[NOTE—The relative retention times for quinidine and dihydroquinidine are 1 and 1.5, respectively.]

Suitability requirements

Resolution: NLT 2.5 between quinidine and dihydroquinidine

Relative standard deviation: NMT 2.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Calculate the percentage of the labeled amount of quinidine sulfate (quinidine sulfate and dihydroquinidine sulfate) in the portion of the Tablets taken:

$$\text{Result} = (r_{B,U} + r_{D,U}) / (r_{B,S} + r_{D,S}) \times (C_S / C_U) \times 100$$

$r_{B,U}$ = peak response of quinidine from the *Sample solution*

$r_{D,U}$ = peak response of dihydroquinidine from the *Sample solution*

$r_{B,S}$ = peak response of quinidine from the *Standard solution*

$r_{D,S}$ = peak response of dihydroquinidine from the *Standard solution*

C_S = concentration of [USP Quinidine Sulfate RS](#) in the *Standard solution* (mg/mL)

C_U = nominal concentration of quinidine sulfate in the *Sample solution* (mg/mL)

Acceptance criteria: 90.0%–110.0%

PERFORMANCE TESTS

• [DISSOLUTION \(711\)](#)

Medium: 0.01 N hydrochloric acid; 900 mL

Apparatus 1: 100 rpm

Time: 30 min

Standard solution: [USP Quinidine Sulfate RS](#) in *Medium*

Sample solution: Dilute with *Medium* to a concentration that is similar to the *Standard solution*. Pass a portion of the solution under test through a suitable filter.

Instrumental conditions

Mode: UV

Analytical wavelength: 248 nm

Analysis

Samples: *Standard solution* and *Sample solution*

Tolerances: NLT 85% (Q) of the labeled amount of quinidine sulfate $[(C_{20}H_{24}N_2O_2)_2 \cdot H_2SO_4 \cdot 2H_2O]$ is dissolved.

Change to read:

- [UNIFORMITY OF DOSAGE UNITS \(905\)](#): ▲ Meet the requirements ▲ (CN 1-Aug-2023)

Procedure for content uniformity

Standard solution: 40 µg/mL of [USP Quinidine Sulfate RS](#) in dilute hydrochloric acid (1 in 100)

Sample solution: Transfer a crushed Tablet to a 250-mL volumetric flask, add 175 mL of dilute hydrochloric acid (1 in 100), and shake by mechanical means for 30 min. Add dilute hydrochloric acid (1 in 100) to volume. Filter a portion of the mixture, discarding the first 20 mL of the filtrate. Dilute quantitatively, if necessary.

Instrumental conditions

Mode: UV

Cell: 1 cm

Analytical wavelength: 345 nm

Blank: Water

Analysis

Samples: *Standard solution*, *Sample solution*, and *Blank*

Calculate the percentage of the labeled amount of quinidine sulfate and dihydroquinidine sulfate, calculated as quinidine sulfate $[(C_{20}H_{24}N_2O_2)_2 \cdot H_2SO_4 \cdot 2H_2O]$, in the Tablet taken:

$$\text{Result} = (A_U / A_S) \times (C_S / C_U) \times 100$$

A_U = absorbance of the *Sample solution*

A_S = absorbance of the *Standard solution*

C_S = concentration of [USP Quinidine Sulfate RS](#) in the *Standard solution* (mg/mL)

C_U = nominal concentration of quinidine sulfate in the *Sample solution* (mg/mL)