


Status: Currently Official on 16-Feb-2025
 Official Date: Official as of 01-Aug-2022
 Document Type: USP Monographs
 DocId: GUID-9E0D287C-FC88-4E05-B30E-EFA487D01CD2_5_en-US
 DOI: https://doi.org/10.31003/USPNF_M79120_05_01
 DOI Ref: nph4c

© 2025 USPC
 Do not distribute

Sulfamethazine



$C_{12}H_{14}N_4O_2S$ 278.33

Benzenesulfonamide, 4-amino-*N*-(4,6-dimethyl-2-pyrimidinyl)-;

*N*¹-(4,6-Dimethyl-2-pyrimidinyl)sulfanilamide CAS RN®: 57-68-1; UNII: 48U51W007F.

DEFINITION

Sulfamethazine contains NLT 99.0% and NMT 100.5% of sulfamethazine ($C_{12}H_{14}N_4O_2S$), calculated on the dried basis.

IDENTIFICATION

• **A. SPECTROSCOPIC IDENTIFICATION TESTS (197), Infrared Spectroscopy:** 197K

• **B. PROCEDURE**

Sample: 0.10 g of Sulfamethazine

Analysis: To the *Sample* add 10 mL of [water](#) and just sufficient 0.1 N [sodium hydroxide](#) (4 mg/mL) to give a faint pink spot on phenolphthalein paper. Add 5 drops of [cupric sulfate TS](#).

Acceptance criteria: A yellow-green precipitate is formed, and it becomes brown on standing.

ASSAY

• **PROCEDURE**

Analysis: Proceed as directed in [Nitrite Titration \(451\)](#).

Calculate the percentage of sulfamethazine ($C_{12}H_{14}N_4O_2S$) in the portion of Sulfamethazine taken:

$$\text{Result} = V \times M \times F \times (1/W) \times 100$$

V = volume of *Titrant* consumed by the *Sample* (mL)

M = actual *Titrant* concentration, in molarity (mEq/mL)

F = equivalency factor of sulfamethazine, 278.3 mg/mEq

W = weight of the *Sample* taken (mg)

Acceptance criteria: 99.0%–100.5% on the dried basis

IMPURITIES

• **RESIDUE ON IGNITION (281):** NMT 0.1%

Delete the following:

▲ **SELENIUM (291):**

Sample: 200 mg of Sulfamethazine

Acceptance criteria: NMT 30 ppm▲ (USP 1-Aug-2022)

• **ORDINARY IMPURITIES (466).**

Sample solution: [Acetone](#)

Standard solution: [Acetone](#)

Eluant: [Ethyl acetate](#), [methanol](#), and [ammonium hydroxide](#) (17:6:5)

Visualization: 11**SPECIFIC TESTS**

- **MELTING RANGE OR TEMPERATURE** (741): 197°–200°

- **ACIDITY**

Sample: 3.0 g of Sulfamethazine**Sample solution:** Digest the *Sample* with 150 mL of [carbon dioxide-free water](#) at 70° for about 5 min, stirring occasionally to maintain suspension of the *Sample*. Cool the mixture rapidly in an ice bath to 20 ± 0.5°, with mechanical stirring. Filter immediately, with vacuum, omitting washing of the solid but drying it thoroughly by suction. Use the filtrate.**Titrimetric system**(See [Titrimetry](#) (541).)**Mode:** Direct titration**Titrant 1:** [0.1 N sodium hydroxide VS](#)**Titrant 2:** [0.1 M sodium nitrite VS](#)**Endpoint detection:** Visual**Analysis 1:** To 25.0 mL of the *Sample solution* add 2 drops of [thymolphthalein TS](#). Titrate with *Titrant 1*. Record the volume of *Titrant 1* consumed (V_1).**Analysis 2:** To a second 25.0-mL portion of the *Sample solution* add 10 mL of [hydrochloric acid](#). Cool the mixture to 15°, and proceed as directed in [Nitrite Titration](#) (451). Record the volume of *Titrant 2* consumed (V_2).

$$\text{Result} = V_1 - V_2$$

 V_1 = volume of *Titrant 1* consumed by the *Sample solution* in *Analysis 1* (mL) V_2 = volume of *Titrant 2* consumed by the *Sample solution* in *Analysis 2* (mL)**Acceptance criteria:** NMT 0.5 mL

- **LOSS ON DRYING** (731).

Sample: 1 g of Sulfamethazine**Analysis:** Dry at 105° for 2 h.**Acceptance criteria:** NMT 0.5%

- **CLARITY AND COLOR OF SOLUTION**

Sample: 1 g of Sulfamethazine**Analysis:** Dissolve the *Sample* in a mixture of 20 mL of [water](#) and 5 mL of 1 N [sodium hydroxide](#).**Acceptance criteria:** The solution is clear and not more deeply colored than pale yellow.**ADDITIONAL REQUIREMENTS**

- **PACKAGING AND STORAGE:** Preserve in well-closed, light-resistant containers.

- **USP REFERENCE STANDARDS** (11).

[USP Sulfamethazine RS](#)**Auxiliary Information** - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
SULFAMETHAZINE	Documentary Standards Support	SM32020 Small Molecules 3
REFERENCE STANDARD SUPPORT	RS Technical Services RSTECH@usp.org	SM32020 Small Molecules 3

Chromatographic Database Information: [Chromatographic Database](#)**Most Recently Appeared In:**

Pharmacopeial Forum: Volume No. 47(2)

Current DocID: GUID-9E0D287C-FC88-4E05-B30E-EFA487D01CD2_5_en-US**DOI:** https://doi.org/10.31003/USPNF_M79120_05_01